

Treatment Guide

MS-65J1



Product positioning: 4-in-1 comprehensive beauty machine that adopts technologies of Ultrasonic nutrient introduction, EMS, Hot&Cold and Eye BIO and can replace all facial treatment in beauty salons.

Packing List

- 1x Ultrasonic facial handle
- 1x Quadrupole photon ems for face lifting
- 1x Photon cold hammer for skin rejuvenation
- 1x Eye bio handle

TechnicalParameters

Power Supply Input:110V-220V 50HZ/60HZ

Photon Ultrasonic Handle:

Frequency:1MHz

Power: 10W

4-Polar Photon BIO Handle:

Frequency:4kHz

Power: 3W

Cold & Hot Hammer:

Power (cold) : 21W

Power (hot) : 18W

COLD:32-39.2°F

HOT: 113°F

Eye EMS Handle:

Frequency:62.5Hz

Power: 4W

Product Size: 13.1 x 10 x 8.8 in. (33.5 x 25.5 x 22.5 cm)

Package Size: 38 x 26 x 19 cm

Net Weight: 2.8 kg


Gross Weight: 3.3 kg

Dos and Don'ts


1. The plug with grounding pin must be used before the instrument is used and the power socket of the instrument is properly grounded.
2. Ensure that the voltage of the instrument is suitable. If the voltage of the local power supply is not stable, we suggest that the user add a voltage regulator with matching power between the mains and the instrument.
3. In order to ensure the treatment effect and normal service life of the instrument, please use the specified parts provided or suggested by the original manufacturer.
4. The instrument should not be placed in damp places or near water sources, and should not be directly exposed to sunlight.
5. Do not place the instrument near a strong heat source, as this may affect the life and normal use of the instrument.
6. Before treatment, please remove all metal objects from the body first, so as to avoid unexpected conditions and influence the curative effect.
7. Please do not use this machine on eyes, thyroid, parathyroid, testicles, abdomen of pregnant women, and pacemakers.
8. Please turn off the power switch of the device, and ensure that the total power supply of the device is turned off before the staff leaves after daily use, so as to ensure the safety of the electrical products.
9. Women in pregnancy, lactation or menstruation can not use this machine.
10. Set the energy at lowest level at the beginning and slowly raise it after the customer adapts to it.
11. The treatment area should be kept moisturized during the operation and do not perform the treatment on dry skin.
12. Use machine or train machine operators in strict accordance with instructions in the user manual.
13. Clean the instrument with normal saline after operation to ensure its cleanliness to guarantee its service life.

Control panel

The control panel is as shown below



1. Control buttons



time and its adjusting buttons



energy level and its adjusting buttons



button for choosing a certain treatment



mode button



start/stop button

2. Detailed operation of EMS

press "FUNC" button and choose EMS



time and its adjusting buttons



energy level and its adjusting buttons

press "MODE" button and choose working mode

M1: constant working mode, suitable for those with dry, rough and dark-toned skin.

M2: high-frequency discontinuous working mode, suitable for first-time users.

M3: low-frequency discontinuous working mode, suitable for those with sensitive skin



start/stop button



Unscrew the upper part of the handpiece, take out the glass container inside and put serum into it. Then screw upper part of the handpiece tight. (please note: serum in high thickness is not allowed)

3. Detailed operation of Ultrasound

press "FUNC" button and choose Ultrasound



time and its adjusting buttons



energy level and its adjusting buttons press "MODE" button and

choose working mode

M1: constant working mode, suitable for those with dry, rough and dark-toned skin.

M2: discontinuous working mode, suitable for first-time users or those with sensitive skin



start/stop button

4. Detailed operation of Hot&Cold

press "FUNC" button and choose Hot&Cold



time and its adjusting buttons



energy level and its adjusting buttons

press “MODE”button and choose working mode

M1: COLD mode, suitable to calm&soothe skin, relieve acne&repair skin and minimize pore size.

M2 : HOT mode, suitable to temporarily expand pore size and help with skin’s absorption.



start/stop button

5. Detailed operation of BIO

press “FUNC”button and choose BIO



time and its adjusting buttons




energy level and its adjusting buttons

press “MODE”button and choose working mode

M1: constant working mode, suitable for those with severe skin problems around eyes.

M2: discontinuous working mode, suitable for first-time users or those with sensitive skin around eyes.



The treatment receiver should hold this before treatment.



start/stop button

The science behind the treatment

EMS

EMS is the abbreviation of electronic muscle stimulation, which is known as muscle electrical stimulation technology, and also known as muscle movement technology. The simple explanation is that under normal circumstances, muscles begin to make movements after receiving signals from the brain, while EMS directly conveys the signal to the muscle through current stimulation, and then causes the muscle to make passive movements so as to achieve the effect of exercise.

Effects

1. It can relax muscles, promote local blood circulation, prevent muscle atrophy, and relieve muscle spasm to exercise muscles. It can accurately and quickly help firm skin for treatment receivers.
2. The low-intensity micro electric current can stimulate muscles to contract, making muscles stronger instead of reducing them. So, it can relieve the sag caused by muscle relaxation and make face contour clearer.

Ultrasound

It makes high-frequency ultrasonic wave penetrate into the subcutaneous tissue at a deep layer of skin. Ultrasonic wave has high energy and strong permeability, which can drive facial cells to vibrate accordingly. It can not only drive pore grease out, but also change the volume of cells to improve blood circulation, which can stimulate cell regeneration and make the skin shiny and elastic. Secondly, because of the heat effect caused by its high-frequency vibration, it also has the function of quickly introducing nutrition into the skin, which can make the skin better absorb the nutrients of the skin care products.

Hot&Cold

It is based on thermal expansion and cold contraction. Our skin reduces the metabolic rate at 5 degrees Celsius, slows down the movement of cells and slows down the sensitivity of nerve endings, which inhibits pain and swelling of skin for allergic reactions. Due to thermal expansion and cold contraction, the temperature of 5 degrees Celsius causes the blood vessels to contract, resulting in pore contraction to restrain thick pores and firm the skin. Therefore, it can smooth facial wrinkles and crow's feet.

Physiological effects of 5°C :

1. It can constrict blood vessels, resulting in tightening of pores, which can greatly improve the size of pores.
2. It can strengthen myofibril, make the skin firm and elastic, and smooth the wrinkles and crow's feet.

3. It can constrict blood vessels, and reduce varicose veins and red veins effectively. It is easier to purify and detoxify venous blood in terms of harmful substances, and it has a significant effect on reduction and elimination of spots.
4. It can reduce the metabolic rate, inhibit and slow down the cell activity, reduce the sensitivity of the peripheral nerve to play a calming and soothing effect, and relieve skin swelling and pain, especially for allergic skin.

Physiological response of 40°C treatment :

1. 40°C makes people feel warm and comfortable, which will be able to relax muscles, keep the balance of vegetative nervous system and relieve pain and soreness. It can promote blood circulation, increase the oxygen content in the blood, and open pores and extend blood vessels to make skin absorb nutrients quickly, which is a milder and more comfortable absorbing way.
2. It can increase the metabolic rate of the skin, increase the heat energy content of the tissue, and make skin naturally eliminate the waste and toxins, which has a certain effect on the elimination of bags under the eyes and dark circles.

The alternating use of cold and heat hammer has multiplying effect on skincare, and can activate cells, and improve and purify the skin, which is a good instrument for health and beauty treatment.

BIO

It stimulates the energy contained in the ATPs of cells to restore the normal operation and function of cells after its imitating human biological current enters into muscle cells through skin to provide energy for muscle and body. It can improve cell permeability to make it 30%-40% higher than usual. It can help the body better absorb nutrition, water, oxygen etc., and eliminate wastes, toxins and retained liquids as soon as possible. It can promote the synthesis of muscle fibers to make it 68% higher than usual, keep the skin moisturized, soft, firm and shiny, and increase collagen synthesis by 73% to reduce wrinkles and leave the skin softer, shinier and more elastic. In this way, it can deliver the results of firming the skin around the eyes, removing bags under the eyes, removing dark circles under the eyes, and relieving bloodshot eyes and fine lines.

Skin-tightening&anti-ageing

1. Needed handpiece(s)

EMS, Ultrasound

2. Effects

1. Firm face and lead nutrients into deep skin.

2. Increase skin elasticity and exquisiteness.
3. Moisturize skin and increase the absorption rate of the skin.
4. Accelerate blood circulation and metabolism.
5. Relieve the condition of double chin and tighten skin.
6. Stimulate collagen regeneration and slow down aging.

3. Applicable Range

1. Those whose faces are dark and dull.
2. Those whose skin is flabby and sagging.
3. Those who have fine lines, nasolabial folds, wrinkles at eye corners.
4. Those whose face contours are not clear.
5. Those who have rough skin and large pores.
6. Those who work for a long time in the office and are exposed to ultraviolet rays.
7. Those who have sagging skin, edema or puffiness after giving birth

4. Inapplicable Range

1. Those who have just undergone plastic surgery with prostheses and metal materials, etc., implanted at the operation area.
2. Those who have used injection products recently, such as hyaluronic acid, skin booster or wrinkle remover.
3. Those who have hypertension, heart disease, diabetes, severe hyperthyroidism, malignant tumor, etc.
4. Those who are having allergic reactions, who have severely sensitive skin or metal allergies.
5. Those who have wounds on skin.
6. Those who are of advanced years.
7. Those who have skin diseases and infectious diseases.

5. Treatment steps(75 min)

Something needed before treatment: washing basin, disposable bag, disposable washcloth, makeup remover, face cleanser, toner, massage cream(essential oil), essence, face mask, lotion/face cream.

1. Remove makeup and clean face(remove the residue of makeup, and clean less stubborn dirt in pores)
2. Evaluate skin: let the treatment candidate look her in the mirror(give personalized skincare treatment according to the skin problems such as uneven skin tone, pigmentation, wrinkle.)

3. Apply toner(clean face again and keep skin moisturized)

4. Massage techniques(relax the treatment receiver. Heat up facial skin and promote circulation. Temporarily expand pore size to help remove grease in pores.(8-15 min)

1. Apply massage cream(massage essential oil) evenly to face and soothe face, 3 times.
2. Apply digital pressing technique on the following acupoints (RN-24, DU-26, ST-4, ST-6, LI-20, BL-1, BL-2, EM3, TE-23, the temple, GB-1, ST-1, ST-2), 3 times.
3. Alternately apply lifting technique from chin to ear lobe, mouth corner to Ermen(SJ21), and nose wing to temple, lift eye corner from lower eyelid to temple with both hands, 3 times. Repeat the technique on the other side.
4. Apply lifting technique towards hairline on forehead, 3 times
5. One side, apply plucking technique with both hands alternately to lift face, 1 min. Repeat the technique on the other side.
6. Move in a zigzag motion on forehead with middle and ring fingers, slide to the front side of ears and lift back and forth for 3 times and then slide to the back of ears.

5. Clean the area and apply essence evenly

6. Ultrasound treatment(help remove grease from pores, promote skin's absorption of nutrients, stimulate cell renewal and leave skin full of glow and elasticity)

- (1). Set parameters: time:10-15 min, energy level: 3-7, mode: M1/M2
- (2). Move the handpiece in spiral curve from chin towards temple, 3-5 times. Do the same on the other side.
- (3). Move the handpiece upwards from chin in a line-by-line way, 3-5 times. Do the same on the other side.
- (4). Move the handpiece in spiral curve towards temples on forehead, 3-5 times.
- (5). Move the handpiece towards hairlines on forehead, 3-5 times.

7. Clean face

8. EMS treatment(micro electric currents can stimulate the movements of facial muscles and restore the functions of cells, thus relieving the condition of saggy skin due to muscle relaxation to tightening skin and rebuild facial curve)(10-15 min)

- (1). Set parameters: time:10-15 min, energy level: 3-7, mode: M1/M2(choose proper mode according to individual skin's situation)
- (2). Move the handpiece from chin to earlobe, mouth corner to ear center and nose wing to temple line by line respectively. In each move, keep the handpiece on skin for 1-2 secs. Perform the technique, 5-8 times.
- (3). Use one hand and the handpiece to push from chin to earlobe, mouth corner to ear center and nose wing to temple line by line respectively. In each move, keep the

handpiece on skin for 1-2 secs. Perform the technique, 3-5 times.

(4). Repeat the techniques above on the other side.

(5). Move the handpiece towards hairlines on forehead, 3-5 times.

9. Check the treatment results (skin improvements like brighter skin tone, tightened & lifted skin and moisturized skin)

10. Clean face and apply face mask (10-15 min)

11. Remove makeup, clean face and apply toner, essence, lotion/face cream and so on.

12. Recommend follow-up treatment

Make a tailored treatment plan according to the feedback from the treatment receiver about expected results. For example: take 2-3 treatment sessions a week. After one course of treatment that includes 12 sessions, skin will be lifted with glow. Two courses will make skin smooth and plump, leaving face contour clearer. Three courses will promote skin metabolism and detox processes, help fade pigmentation and skin darkness, slow down skin ageing and prevent skin from sagging, and stimulate collagen regeneration to leave skin smooth and firm with glow.

6. Do's and Don'ts after Treatment

1. Wash face with hot water within 3 days
2. Keep hydrated and prevent skin damages caused by sun ray
3. Do not use irritating products within one week such as AHA, scrubbing and exfoliating products
4. Do not stay up late, smoke cigarettes and drink alcohol. Eat more fruits, vegetables and light foods.
5. Apply face mask at least 3 times in the following week.

Skin whitening & rejuvenation

1. Needed Handpieces

EMS, Ultrasound

2. Results

1. Help with the introduction of nutrients to quickly nourish skin.

2. Whiten and rejuvenate skin and relieve dark skin tone.
3. Stimulate the regeneration of cells and increase skin's elasticity.
4. Improve skin texture and increase skin's resistance.
5. Relieve skin's roughness to smooth skin and alleviate the condition of large pore size.

3. Applicable range

1. Those who have dry skin and whose skin is less nourished.
2. Those with rough skin and large pores.
3. Those with uneven-toned, dark-toned skin.
4. Those who have poor skin absorption and lack moisture at the basal layer of skin.

4. Inapplicable range

1. Women in menstruation, pregnancy or lactation.
2. Those who "3Hs", hypothyroidism, malignant tumors and other severe diseases.
3. Those who have just taken microneedling, skin booster injections, invasive plastic surgery and have unhealed incisions.
4. Those with skin disease or infectious diseases.
5. Those who are suffering from severe allergies(like rashes)

5. Treatment steps(75 min)

Something needed before treatment: washing basin, disposable bag, disposable washcloth, makeup remover, face cleanser, toner, massage cream(essential oil), skin-whitening essence, skin-whitening face mask, lotion/face cream.

1. Remove makeup and clean face(remove the residue of makeup, and clean less stubborn dirt in pores)
2. Evaluate skin: let the treatment candidate look her in the mirror(give personalized skincare treatment according to the skin problems such as uneven skin tone, pigmentation, wrinkle.)
3. Apply toner(clean face again and keep skin moisturized)
4. Massage techniques(relax the treatment receiver. Heat up facial skin and promote circulation. Temporarily expand pore size to help remove grease in pores.(8-15 min)
 1. Apply massage cream(massage essential oil) evenly to face and soothe face, 3 times.
 2. Apply digital pressing technique on the following acupoints (RN-24, DU-26, ST-4,

ST-6, LI-20, BL-1, BL-2, EM3, TE-23, the temple, GB-1, ST-1, ST-2), 3 times.

3. Alternately apply lifting technique from chin to ear lobe, mouth corner to Ermen(SJ21), and nose wing to temple, lift eye corner from lower eyelid to temple with both hands, 3 times. Repeat the technique on the other side.

4. Apply lifting technique towards hairline on forehead, 3 times

5. One side, apply plucking technique with both hands alternately to lift face, 1 min. Repeat the technique on the other side.

6. Move in a zigzag motion on forehead with middle and ring fingers, slide to the front side of ears and lift back and forth for 3 times and then slide to the back of ears.

5. Clean the area and apply essence evenly

6. Ultrasound treatment(help remove grease from pores, promote skin's absorption of nutrients, stimulate cell renewal and leave skin full of glow and elasticity)

(1). Set parameters: time:10-15 min, energy level: 3-7, mode: M1/M2

(2). Move the handpiece in spiral curve from chin towards temple, 3-5 times. Do the same on the other side.

(3). Move the handpiece upwards from chin in a line-by-line way, 3-5 times. Do the same on the other side.

(4). Move the handpiece in spiral curve towards temples on forehead, 3-5 times.

(5). Move the handpiece towards hairlines on forehead, 3-5 times.

7. Clean face

8. EMS treatment(micro electric currents can stimulate the movements of facial muscles and restore the functions of cells, thus relieving the condition of saggy skin due to muscle relaxation to tightening skin and rebuild facial curve)(10-15 min)

(1). Set parameters: time:10-15 min, energy level: 3-7, mode: M1/M2(choose proper mode according to individual skin's situation)

(2). Move the handpiece from chin to earlobe, mouth corner to ear center and nose wing to temple line by line respectively. In each move, keep the handpiece on skin for 1-2 secs. Perform the technique, 5-8 times.

(3). Use one hand and the handpiece to push from chin to earlobe, mouth corner to ear center and nose wing to temple line by line respectively. In each move, keep the handpiece on skin for 1-2 secs. Perform the technique, 3-5 times.

(4). Repeat the techniques above on the other side.

(5). Move the handpiece towards hairlines on forehead, 3-5 times.

9. Check the treatment results(skin improvements like brighter skin tone, tightened&lifted skin and moisturized skin)

10. Clean face and apply face mask(10-15 min)

11. Remove makeup, clean face and apply toner, essence, lotion/face

cream and so on.

12. Recommend follow-up treatment

Make a tailored treatment plan according to the feedback from the treatment receiver about expected results. For example: take 2-3 treatment sessions a week. After one course of treatment that includes 12 sessions, skin will lifted with glow. Two courses will make skin smooth and plump, leaving face contour clearer. Three courses will promote skin metabolism and detox processes, help fade pigmentation and skin darkness, slow down skin ageing and prevent skin from sagging, and stimulate collagen regeneration to leave skin smooth and firm with glow.

6. Do's and Don'ts after Treatment

1. Wash face with hot water within 3 days
2. Keep hydrated and prevent skin damages caused by sun ray
3. Do not use irritating products within one week such as AHA, scrubbing and exfoliating products
4. Do not stay up late, smoke cigarettes and drink alcohol. Eat more fruits, vegetables and light foods.
5. Apply face mask at least 3 times in the following week.

Skin soothing and repairing

1. Needed handpieces

EMS, Hot&Cold handpiece

2. Results

1. Quickly calm skin and efficiently keep skin moisturized.
2. Relieve the condition of sensitive skin.
3. Relieve excess grease on skin and minimize pore size.
4. Soothe skin after photon or micro-current treatment.
5. Relieve acne.

3. Applicable range

1. Those with dry skin.
2. Those who have sensitive skin, and are prone to skin redness, swelling and allergies.
3. Those with large pore size, rough skin, and peeling skin.

4. Inapplicable range

1. Those who have just taken invasive plastic surgery with incisions unhealed.
2. Those with skin disease or infectious diseases.
3. Those who have unrealistic expectations about the results.
4. Those who have just taken injects like hyaluronic acid, skin booster, and wrinkle remover

5. Treatment steps(60 min)

Something needed before treatment: washing basin, disposable bag, disposable washcloth, makeup remover, face cleanser, toner, massage cream(essential oil), essence, face mask, lotion/face cream.

1. Remove makeup and clean face(remove the residue of makeup, and clean less stubborn dirt in pores)
2. Evaluate skin: let the treatment candidate look her in the mirror(give personalized skincare treatment according to the skin problems such as uneven skin tone, pigmentation, wrinkle.)
3. Apply toner(clean face again and keep skin moisturized)
4. Apply restoration essence.
5. Cold&Hot handpiece treatment(calm and soothe skin, relieve acne and repair skin, minimize pore size) (5-8 min)
Set parameters: time: 5-8 min, energy level: 3-7, mode: M1
 - (1). Move the handpiece in spiral curve from chin towards temple, 5-8times. Do the same on the other side.
 - (2). Move the handpiece in spiral curve towards temples on forehead, 5-8 times.
 - (3). Repeat the techniques 3-5 times on the area where the condition of skin redness or acne is severe.
6. EMS treatment(micro electric currents can stimulate the movements of facial muscles and restore the functions of cells, thus relieving the condition of saggy skin due to muscle relaxation to tightening skin and rebuild facial curve)(10-15 min)
 - (1). Set parameters: time:10-15 min, energy level: 3-7, mode: M2/M3(choose proper mode according to individual skin's situation)
 - (2). Move the handpiece from chin to earlobe, mouth corner to ear center and nose wing to temple line by line respectively. In each move, keep the handpiece on skin for 1-2 secs. Perform the technique, 5-8 times.
 - (3). Repeat the techniques above on the other side.
 - (4). Move the handpiece towards hairlines on forehead, 3-5 times.

7. Check the treatment results(skin improvements like brighter skin tone, tightened&lifted skin and moisturized skin)

8. Clean face and apply face mask(10-15 min)

9. Remove makeup, clean face and apply toner, essence, lotion/face cream and so on.

10. Recommend follow-up treatment

Make a tailored treatment plan according to the feedback from the treatment receiver about expected results. For example: take 2-3 treatment sessions a week. After one course of treatment that includes 12 sessions, skin will be return to normal and less prone to allergic reactions caused by seasonal factors, cosmetic products and dryness of skin. After two courses, skin will be supple and velvety. There will be ruddy cheeks thanks to promoted blood circulation. After three courses, effects will be consolidated.

Wrinkle removal for eye area

1. Needed handpiece

BIO, the metal ring

2. Results

1. Relieve skin hyperpigmentation and fade dark circles.
2. Relieve under-eye bags and edema around eyes.
3. Promote circulation around eyes, relieve Dry Eye Syndrome(DES)
4. Promote metabolism and prevent hyperpigmentation
5. Provide eye area with nutrients to nourish skin.
6. Reduce dry lines and fine lines around eyes.

3. Applicable Range

1. Those who have dry skin around eyes and whose skin lacks nutrients around eyes.
2. Those who often stay up late and have dark circles.
3. Those have have dry lines and fine lines around eyes due to dry skin.
4. Those who suffer from under-eye bags and edema around eyes.

5. Those who have poor skin absorption around eyes.

4. Inapplicable range

1. Women in menstruation, pregnancy or lactation.
2. Those who “3Hs”, hypothyroidism, malignant tumors and other severe diseases.
3. Those who have just taken microneedling, skin booster injections, invasive plastic surgery and have unhealed incisions.
4. Those with skin disease or infectious diseases.
5. Those who are suffering from severe allergies(like rashes)

5. Treatment steps(45min)

Something needed before treatment: washing basin, disposable bag, disposable washcloth, makeup remover, face cleanser, toner, eye essence, eye mask, eye cream.

1. Remove makeup and clean face(remove the residue of makeup, and clean less stubborn dirt in pores)
2. Evaluate skin: let the treatment candidate look her in the mirror(give personalized skincare treatment according to the skin problems such as uneven skin tone, pigmentation, wrinkle.)
3. Apply toner(clean face again and keep skin moisturized)
4. Massage techniques(relax the treatment receiver. Heat up skin around eyes and promote circulation.(5-8 min)
 - (1). Apply eye essence evenly, and soothe eye area 3 times.
 - (2). Apply digital pressing technique on the following acupoints(BL-1, BL-2, EM3, TE-23, the temple, GB-1, ST-1, ST-2), 3 times.
 - (3). Use middle finger to lift from inner eye corner to temple back and forth, 3-5 times. Repeat the technique on the other side.
 - (4). Use middle finger and ring finger to lift from inner eye corner to temple along eyebrow, 3-5 times. Repeat the technique on the other side.
5. Clean the area and apply eye essence evenly
6. BIO treatment(relieve under-eye bags, dark circles, fine lines and the condition of bloodshot eyes) (10 min)

Let the treatment receiver hold the metal ring.
Set parameters: time: 10 min, energy level: 3-7, mode: M1

 1. Move the handpiece to lift skin from lower eyelid to temple in small circles, 3-5 times.
 2. Move the handpiece to lift skin from lower eyelid to eye corner, 3-5 times.

3. Move the handpiece to lift skin from lower eyelid to temple, 3-5 times.
4. Use one hand and the handpiece to push from lower eyelid to temple. 3-5 times.
5. Move the handpiece to lift skin from eyebrow towards hairlines, 3-5 times.
6. Repeat the techniques on the other side
7. Check the treatment results.
8. Clean the area and apply eye mask(10-15 min)
9. Remove eye mask, clean the area and apply toner, eye essence, eye cream.
10. Recommend follow-up treatment

Make a tailored treatment plan according to the feedback from the treatment receiver about expected results. For example: take 2-3 treatment sessions a week. After one course of treatment that includes 12 sessions, fine wrinkles and black circles are reduced and skin colour is brightened. After three courses, skin in eye area is totally firmed and rejuvenated with glossiness. If the customers keep using this machine, eye aging can be relieved and prevented.