

# User Guide

## MS-12Y1



# Preface

Dear users:

Thanks for choosing our company's product, in order to fully make use of the products, we sincerely suggest you to follow:

Please read the instruction in detail and keep it well for further reading and consulting.

Please follow the instruction to fit and operate the machine correctly. Please do not remove and change any accessories of the machine. It forbids others to open and remove the mainframe except the authorized people in the company. It forbids others to change and maintain the machine except the professional people. Please contact us if there is any trouble during operation.

After Sale service

We offer ONE YEAR warranty. Within warranty period, all exchanges or replacement due to

defects are free. (Buyer is responsible for shipping charges).

If the machine is over warranty, please also do not hesitate to contact us. we will assure you of our best service all the time.

This product does not belong to medical, please don't use for disease treatment!

# Table of Contents

Preface
Part I
Brief Introduction
Inapplicable Range
Part II
1. Detailed Operations
2. Technical Parameters
3. Safety Precautions
4. Dos and Don'ts
5. Potential Side Effects
6. Packing List
7. Operation Steps

# Part I

## Brief Introduction

RF system beauty system uses an advanced new RF Radio Frequency facial and body technology, which can reach the efficacy as following:

**Managing wrinkles or skin rejuvenation:** This device transfer high frequency RF energy to skin and stimulates to produce collagen in skin layer. In this case, stimulated collagen tissues are activated and thus helps prevention of wrinkles due to stress and help expedition of skin rejuvenation.

**Managing freckle and skin whitening:** This device could alleviate freckles due to deteriorated blood circulation in a way that high frequency RF energy is transferred to skin and blood circulation is expedited and it also help skin whitening effects by reducing melanic pigment secretion.



## Inapplicable Range

1. Those with open wounds or lesions on the face and/or neck.
2. Those with severe or cystic acne on the face and/or neck
3. Those with mechanical implants, dermal fillers, implanted electrical devices and metal stents in the face and/or neck area.

# Part II

## 1. Detailed Operation

The detailed interface is as follows:





Quadrupole 3D SMART RF Handle with Photon



3-Polar 3D SMART RF Handle with Photon

Intensity Indicator light: 1-5 Level

Mode: There are two working mode: continuous and Intermittent

Continuous means the energy keep on the same level when using

Intermittent means energy will change from low to high you can change working mode via click "mode" button

ON/OFF: Power on/off , it is connected with the power supply to be charged

## 2. Technical Parameters

Power Supply:110V-220V± 10% 50/60Hz, 12A

Power : 35W

3-Polar 3D Smart RF Handle With Photon: 3Mhz

Quadrupole 3D Smart RF Handle With Photon: 3Mhz

### 3. Safety Precautions

1. Treatment energy is not recommended for use directly on an existing keloid.
2. The RF system has not been evaluated for use in patients on an anticoagulant treatment plan.
3. It is recommended that the following areas should be avoided during treatment:
  - Thyroid gland, thyroid cartilage and trachea
  - Major vessels
4. The RF system has not been evaluated for use in the following patient populations:
  - Pregnant or breast feeding women
  - Children
5. Those with the following disease states:
  - A hemorrhagic disorder or hemostatic dysfunction
  - An active systemic or local skin disease that may alter wound healing
  - Herpes simplex
  - Autoimmune disease
  - Diabetes
  - Bell's palsy

### 4. Dos and Dont's

1. RF system should not be used on a patient's eyes or in a location or technique where ultrasound energy can reach the eye
2. Use this system only if you are trained and qualified to do so.
3. If any problems occur during system operation, ask immediate action(s): lift the transducer off the patient's skin, press the STOP button to completely halt system operation.

### 5. Potential Side Effects

Side effects reported in the clinical evaluation of the RF system were mild and transient in nature. These were limited to:

1. Erythema (redness): The treated area may exhibit erythema immediately following treatment. This typically resolves within a few hours of treatment.
2. Edema (swelling)rythema (redness): The treated area may exhibit erythema immediately following Treatment. This typically resolves within a few hours of treatment.
3. Pain: Momentary discomfort may be experienced during the procedure while energy is being deposited. Post procedure discomfort or tenderness to the touch is also possible.

4. Bruising: Mild bruising, which is caused by damage to soft tissue blood vessels, may occur occasionally and typically resolves within a few days of treatment.
5. Nerve Effects: Transient local muscle weakness may result after treatment due to inflammation of a motor nerve. Transient numbness may result after treatment due to inflammation of a sensory nerve. Transient pain, paresthesia and/or tingling may be experienced No permanent injuries to facial nerves have been reported
6. Scarring: The possibility for scar formation ( which may respond to medical care) May exist if incorrect treatment technique is used

## 6.Packing List

- 1 x Main Machine
- 1 x 3-Polar 3D Smart RF Handle With Photon
- 1x Quadrupole 3D Smart RF Handle With Photon
- 1 x Power Supply Cord

## 7.Operation Steps

Parameter Adjustment	Product	Steps	Diagrams
<b>Tighten Sagging Skin: 30 min, 3 times/week</b>			
Choose function & levels according to your skin type	Makeup remover+ Face wash+ Toner+ Device	<ol style="list-style-type: none"> <li>1.Deep clean the skin Remove make up, wash face, facial tone and then do facial spray to open skin pores.</li> <li>2. Connect the power adapter to main machine</li> <li>3. Plug cable to main machine and connect the other end to the handle</li> <li>4. Press on/off button to start work .press “intensity” button to choose energy level. Click “mode” button to choose work mode</li> <li>5. Put some gel in treatment</li> </ol>	No diagrams

		area before treatment 6. Move the handle and about few minutes later you will feel warm 7. Press ON/Off button again to stop working remove handle and press red on off switch.	
<p><b>Tips:</b> Please keep the three polar RF probe in full contact with the skin during treatment, move the probe on the surface of the cared area, in an even and thorough way.</p> <p><b>Important:</b> The most important thing is to choose suitable cosmetics for your skin. It is not only for therapy, but also an important foundation to provide nutrient to your skin. Therefore, be careful to select the right cosmetics.</p>			