



MS-22Y1

Instructions

Warning

Dear users:

Thanks for choosing our company's product, in order to fully make use of the products, we sincerely suggest follow:

- ◆◆Please read the instruction in detail and keep it well for further reading and consulting.
- ◆◆Please follow the instruction to fit and operate the machine correctly.
- ◆◆Please do not remove and change any accessories of the machine.
- ◆◆It forbids others to open and remove the mainframe except the authorized people in the company.
- ◆◆It forbids others to change and maintain the machine except the professional people.
- ◆◆Please contact us if there is any trouble during operation.

After Sale service

◆We offer ONE YEAR warranty. Within warranty period, all exchanges or replacement due to defects are free. (Buyer is responsible for shipping charges).

If the machine is over warranty, please also do not hesitate to contact us. we will assure you of our best service all the time.

◆How to contact us?

Phone:(+86) 189 -487-55929

Tel: 0755-86368500

Email: sales@mychway.com

This product does not belong to medical, please don't use for disease treatment!

Chapter 1: Safety Warnings

I. Electricity & Machinery Safety

- 1.The machine is either 110V And 220V. So whatever the voltage u can use before starting the machine.
2. Please keep the voltage at stable condition.

II. Contraindication Safety

1. People with heart disease, high blood pressure, or configured cardiac Pacemaker.
2. People with acute inflammation, asthma, deep vein thrombosis, thyroncus, cancer.
3. People with hemorrhagic diseases, trauma or who is bleeding.
4. Pregnant women
5. Medical plastic artificial or metal parts inside the body such as metal teeth or silicone breast implants.
6. People with abnormal immune system.
7. People with numb or insensitive to heat.
8. If you are feeling unwell, Prohibited in the ears, nose, eye throat and do not use this device in advance and tell your doctor.
9. Children under 12 years old are not allowed.

III. Packing Including:

Diamond Dermabrasion Machine

PVC vacuum hose 2pcs

Diamond wands 3pcs

Diamond tips 9 pcs

Cotton filter 1pack

Spray Gun 1pc

Power Adapter 1pc

Chapter 2: How to install accessories.

1. Install filter and rubber band to wands.

Please choose one diamond tip, fix one piece rubber band on the tube, put one piece cotton filter and install the diamond tip
Please connect the tube with the hose



2. Install hose to machine hole and connect hose, please noticed that the outlet for spray gun, the inhale diamond tube.



3. Connect spray gum with machine, please install it like below:



3. After install the accessories, please put them well on the holder

5. Connect power plug with the machine.



6. Press the power button to start/finish treatment.








Chapter 3: program setting

I. Turn the silver knob to adjust the suction of vacuum. Turn to right of the machine, make the vacuum suction at the highest level, turn to left of the machine, vacuum suction is the lowest level.



II. Choose the suitable microdermabrasion tips according to different treatment area and start the treatment, massage along the skin texture.

There are 9 different microdermabrasion tips for you choose:

1		Body, Especially for acne, knee, elbow areas and striae of pregnancy.
2		Face. For the large blemishes.
3		For Scars, fine lines and imperfections.
4		Used on the eyes and both sides of nose. For wrinkles, oil tablets, blackheads and acne.
5		After the D tips peeling procedure, use this smooth tip to vacuum/clear away the dead abraded skin.

III. Operation guide

Before use, clean face with lotion and finish massage face to accelerate blood circulation.



Operating Guide



STEP 1

Cleaning the facial skin and drying it.



STEP 2

Starting from the lower jaw, moving the rejuvenation from bottom to top.



STEP 3

Placing the rejuvenation instrument at the cheeks and moving from the inside to the outside.



STEP 4

Placing the rejuvenation instrument in the forehead area, moving from the middle to both sides.



STEP 5

Placing the rejuvenation instrument in the T zone, moving from top to bottom.



STEP 6

Using Toner to moisturize the skin, skin care products absorbing better.

IV. Use the spary gum to make skin moisturize after use the diamond dermabrasion

1. Use your finger press like below



2. Spray down from customer's forehead. Let customer close eyes and mouth avoiding liquid into, at the same time do not let liquid flow into customer's nose. Please keep the spray gun from face in proper distance.

Chapter 4: After using

- I.** Do a hydrating mask after the treatment.
- II.** Next treatment should be after half of month.