


User Manual

Home SPA Instrument




Content

Packing List.....	1
Technical Specification.....	1
Precautions.....	3
Contraindications.....	4
Interface Introduction.....	5
CAVI 30K Handle Test.....	15
CAVI 50K Handle Test.....	16
CAVI 80K Handle Test.....	16
VACUUM&RF Handle Test.....	17
LIPO LASER Handle Test.....	17
Operation Guide.....	18
Contact Us.....	36

Packing List



Main unit*1




CAVI 30K/50K&
RF handle*1




CAVI 80K handle*1



VACUUM&RF handle*1



Laser pad*6



Handle holder*1 set



Power cord * 1



Fuse*2





Strap*1 set






Filter*1

Technical Specification

	Input Voltage: 110V-260V 50-60HZ		Net Weight: 16.5KG	
	Power: 70W		Gross Weight: 28.5KG	
	Dimensions: 21.7×17.9×16.5inches			
	Input Power: 70W			

	RF Frequency: 1MHz		RF Frequency: 1MHz	
	30K Frequency: 30KHz		50K Frequency: 50KHz	
	Power: 70W		Power: 60W	

	Frequency: 80KHz			Wavelength: 65 0 nm160mw	
	Power: 33W			Power: 31W	

	RF Frequency: 1MHz		Negative pressure flow: 10LPM	
	Power: 68W		50K Frequency: 50KHz	
	Suction intensity: -80 KPa		Noise level (30cm away): <60db (30cm away)	

Precautions

- 1.This product is guaranteed by the designated after-sales service unit.
- 2.Before using the instrument, use a plug with a grounding pin and ensure the power socket is grounded.
- 3.Clean the instrument with normal saline for hygienic reasons and to prolong its service life.
- 4.Always use the designated accessory provided or recommended by the manufacturer to ensure performance throughout its service life.
- 5.This instrument is a high-precision electronic device that requires proper storage. Please avoid storing it in humid conditions, direct sunlight, or near flammables and explosives to prevent damage.
- 6.Do not place the instrument in high-temperature conditions, as this may negatively impact its service life and performance.
- 7.Do not apply the instrument on eyes, thyroid, parathyroid, testicles, pregnant women's abdomen, pacemakers, etc.
- 8.Remove all metal objects from the user's body to prevent unpredictable conditions or reduce treatment effectiveness.
- 9.Please power off the instrument when not in use, and make sure to turn off the general power of the instrument before leaving after daily use to ensure the safety of the instrument.
- 10.Individuals with a metal bracket (excluding IUD) or breast implants installed should avoid using this instrument.
- 11.Individuals who have had a cesarean section or surgery within the last six months, a miscarriage within the last three months, or a natural birth within the last two months should

refrain from using this instrument.

12.Ensure the filter cotton is installed inside the vacuum handle, as the foreign body may block the handle.

13.Always begin the treatment at the lowest energy level and gradually increase it as the user becomes accustomed to the sensation.

14.The use of a medium on the treatment area is necessary, and performing the operation on the treatment area without the medium should be avoided.

15.Do not use the 30K/50K/80K handle on the head, breast, thorax, heart, and back.

16.In case you are taking other anti-adiposity drugs, it is recommended to stop the drugs for 1 to 2 months before treatment. If you want to take the treatment immediately, the course of treatment should be extended.

17.Strictly follow the instructions in the manual to use the instrument or train the operator.

Contraindications

1.Individuals diagnosed with diabetes, heart disease, a history of scarring, or hyperthyroidism.

2.Those who have recently undergone surgery or are in the recovery phase.

3.Individuals with malignant tumors, hemophilia, or hemorrhagic conditions.

4.Individuals with pacemakers, vascular stents, or metal implants, as well as those using heart rate devices or other active electronic devices.

5.Individuals afflicted with skin disease or infectious conditions.

6.Individuals experiencing skin irritation, sensitive skin, or metal allergies.


7.Individuals who have recently undergone plastic surgery or have implants in the treatment area.

8.Individuals who have recently received injections such as Hyaluronic Acid and botulinum, undergone chemotherapy, or plastic surgery should refrain from using this instrument on their face and neck.


9.Menstruating, pregnant, or breastfeeding individuals.

Interface Introduction


Connect the power source, power on the instrument, and the Welcome Interface will subsequently appear.




1.AUTO


Tap  to access the following interface.






 Back to the previous step


1.1.ARMS SHAPING


Tap  to access the following interface.






 Process of the treatment


 Main menu  Reset button  Press the on/off button to initiate the CAVI 80K treatment. The parameters are preset to default settings, and you can adjust the time and energy level based on individual tolerance.


 80K mode (The default mode is M1, and you may switch to another mode manually)
M1: Continuous working mode
M2: Intermittent working mode






 Process of the treatment


 Main menu  Reset button  Press the on/off button to initiate the VACUUM&RF treatment. The parameters are preset to default settings, and you can adjust the time and energy level based on individual tolerance.


 VACUUM mode (The default mode is M1, and you may switch to another mode manually)
M1: Only suction.
M2: At a slow intermittent frequency
M3: At a moderate intermittent frequency
M4: At a fast intermittent frequency
M5: At the fastest intermittent frequency



 Process of the treatment




 Main menu  Reset button  Press the on/off button to initiate the LIPO LASER treatment. The parameters are preset to default settings, and you can adjust the time and energy level based on individual tolerance.


 LIPO LASER mode (The default mode is M1, and you may switch to another mode manually)
M1: The lights are steadily on.
M2: The lights flash slowly.
M3: The lights flash quickly.



The entire ARMS SHAPING treatment is considered complete when this concluding interface appears.

1.2.WAIST SHAPING


Tap    to access the following interface.



Process of the treatment


 Main menu

 Reset button

 Press the on/off button to initiate the CAVI 30K&RF treatment. The parameters are preset to default settings, and you can adjust the time and energy level based on individual tolerance.

MODE  M2

30K mode (The default mode is M1, and you may switch to another mode manually)
M1: Continuous working mode
M2: Intermittent working mode




Process of the treatment

 Main menu

 Reset button



Press the on/off button to initiate the VACUUM&RF treatment. The parameters are preset to default settings, and you can adjust the time and energy level based on individual tolerance.

VACUUM MODE  M2 M3 M4 M5

VACUUM mode (The default mode is M1, and you may switch to another mode manually)


M1: Only suction.

M2: At a slow intermittent frequency

M3: At a moderate intermittent frequency

M4: At a fast intermittent frequency

M5: At the fastest intermittent frequency



Process of the treatment

 Main menu

 Reset button



Press the on/off button to initiate the LIPO LASER treatment. The parameters are preset to default settings, and you can adjust the time and energy level based on individual tolerance.


MODE  M2 M3

LIPO LASER mode (The default mode is M1, and you may switch to another mode manually)

M1: The lights are steadily on.


M2: The lights flash slowly.


M3: The lights flash quickly.



The entire WAIST SHAPING treatment is considered complete when this concluding interface appears.

1.3.LEGS SHAPING

Tap  to access the following interface.



Process of the treatment



Main menu



Reset button



Press the on/off button to initiate the CAVI 50K&RF treatment. The parameters are preset to default settings, and you can adjust the time and energy level based on individual tolerance.



50K mode (The default mode is M1, and you may switch to another mode manually)
M1: Continuous working mode
M2: Intermittent working mode



Process of the treatment




Main menu



Reset button



Press the on/off button to initiate the VACUUM&RF treatment. The parameters are preset to default settings, and you can adjust the time and energy level based on individual tolerance.



VACUUM mode (The default mode is M1, and you may switch to another mode manually)


M1: Only suction.

M2: At a slow intermittent frequency

M3: At a moderate intermittent frequency

M4: At a fast intermittent frequency

M5: At the fastest intermittent frequency



Process of the treatment



Main menu



Reset button



Press the on/off button to initiate the LIPO LASER treatment. The parameters are preset to default settings, and you can adjust the time and energy level based on individual tolerance.



LIPO LASER mode (The default mode is M1, and you may switch to another mode manually)


M1: The lights are steadily on.

M2: The lights flash slowly.

M3: The lights flash quickly.




Tap this button to access the LIPO LASER operation interface early, and tap the on/off button to start the LIPO LASER treatment.




The entire WAIST SHAPING treatment is considered complete when this concluding interface appears.


2.MANUAL








Tap  to access the following interface.




2.1.CAVI 30K&RF


Tap  to access the following interface.










-  Time display and set buttons (0 - 60 minutes, with 15 minutes recommended)
-  CAVI 30K energy display and set buttons (range 0-10, 3-8 recommended)
-  RF energy display and set buttons (0-10 levels, with 3-8 being recommended)
-  30K working mode M1: Continuous working mode M2: Intermittent working mode
-  Main menu  Reset button  On/off button


2.2.CAVI 50K&RF


Tap  to access the following interface.



-  Time display and set buttons (0 - 60 minutes, with 15 minutes recommended)
-  CAVI energy display and set buttons (range 0-10, 5-8 recommended)
-  RF energy display and set buttons (0-10 levels, with 3-8 being recommended)
-  50K working mode M1: Continuous working mode M2: Intermittent working mode
-  Main menu  Reset button  On/off button

2.3.VACUUM&RF

Tap  to access the following interface.



< 10 : 00 >
TIME

< 10 >
RF

< 2.0 >
VACUUM SUCTION

< 0.8 >
VACUUM RELEASE

VACUUM MODE

M1 M2 M3 M4 M5

Vacuum mode
M1: Only suction. M2: At a slow intermittent frequency
M3: At a moderate intermittent frequency
M4: At a fast intermittent frequency
M5: At the fastest intermittent frequency


Main menu


Reset button

On/off button

Suction control knob (set according to individual tolerance)

2.4.CAVI 80K

Tap  to access the following interface.



< 10 : 00 >
TIME

< 10 >
CAVI 80K

< 8 >
RF

MODE

M1 M2


80K working mode M1: Continuous working mode M2: Intermittent working mode


Main menu

Reset button

On/off button

2.5.LIPO LASER

Tap  to access the following interface.



Display of the LIPO LASER energy

Energy level set buttons (0 - 10 levels, with 3 - 8 being recommended)

< 20 : 00 >
TIME

Time display and set buttons (0 - 60 minutes, with 15 minutes recommended)

MODE

M1 M2 M3

LIPO LASER working mode M1: Lights are steadily on.
M2: Lights flash slowly. M3: Lights flash quickly.

Main menu

Reset button

On/off button

CAVI 30K HANDLE TEST

- 1.Power on the instrument.
- 2.Select "CAVI 30K&RF" on the main menu.
- 3.Set the energy level to 10.
- 4.Tap the on/off button and drop some water on the treatment head.
- 5.Inspect the water on the treatment head. If it is atomized, the CAVI 30K is functioning normally. If not, please contact the dealer for assistance.

CAVI 50K HANDLE TEST

- 1.Power on the instrument.
- 2.Select "CAVI 50K&RF" on the main menu.
- 3.Set the energy level to 10.
- 4.Inspect the water on the treatment head. If it is vibrated, the CAVI 50K is functioning normally. If not, please contact the dealer for assistance.
- 5.Inspect the water on the treatment head. If it is vibrated, the CAVI 50K is functioning normally. If not, please contact the dealer for assistance.

Note: Please note that the water appearance differs during the 30K and 50K tests. The water is well atomized in the 30K test, whereas visible ripples appear in the water during the 50K test.

CAVI 80K HANDLE TEST

- 1.Power on the instrument.
- 2.Select "CAVI 80K" on the main menu.
- 3.Set the energy level to 10.
- 4.Tap the on/off button and drop some water on the treatment head.
- 5.Inspect the water on the treatment head. If there are visible ripples in the water, the 80K is functioning normally. If not, please contact the dealer for assistance.

VACUUM&RF HANDLE TEST

- 1.Power on the instrument.
- 2.Select " VACUUM&RF" on the main menu.
- 3.Adjust the mode to M1, and set the suction intensity to the maximum using the suction control knob.
- 4.Tap the on/off button and place the handle on your palm.
- 5.If the handle successfully suctions your palm, the VACUUM&RF is operating normally. If not, please contact the dealer for assistance.

LIPO LASER HANDLE TEST

- 1.Power on the instrument.
- 2.Select "LIPO LASER" on the main menu.
- 3.Set the energy level to 10.
- 4.Tap the on/off button and verify if the lights illuminate. If they do, then the LIPO LASER is operating normally. If not, please contact the dealer for assistance.

Operation Guide

Arm Refinement

Handles: CAVI 80K, VACUUM&RF, and Laser Pad

Duration: 60 minutes

Preparation: hot towel, essential oil, gel, strap.

1.Gently wipe the arms with a hot towel.

2.Manual massage (clear channels, enhance circulation)

Product required: essential oil

Recommended working time: 10 minutes

Recommended application:

a.Apply essential oil evenly to the treatment area.

b.Place the arms flat, and use the first web of both hands to push the outer side of the arm.

c.Soothing act.

d.Warm the outer side of the arm by rubbing it with bent fingers.

e.Soothing act.

f.Position the inner side of the arm upward, and push it towards the armpit using the first web of both your hands.

g.Soothing act.

h.Warm the inner side of the arm by rubbing it with bent fingers.

i.Soothing act.

j.Repeat the same on the other arm.

k.Wipe the treatment area with a towel.

3.CAVI 80K - clear channels, enhance body contours

Product required: Gel

Parameters setting:

Recommended working time: 15 minutes

Recommended energy: 3

Recommended Mode: M1

Recommended application:

a.Apply gel evenly to the treatment area.

b.Lift slowly from the lower arm towards the elbow, and from the elbow towards the shoulder using the handle.

c.Slowly move the handle in circular motions from the elbow towards the shoulder.

d.Slowly lift the upper arm from the inside outward.

e.Move the handle in circles on the upper arm.

f.Repeat on the other side.

g.Wipe arms with a towel.

4.VACUUM&RF - Clarify channels, tighten skin, and sculpt the body

Product required: essential oil

Parameters setting:

Recommended working time: 15 minutes

Recommended RF energy: 3

Recommended suction duration: 0.5

Recommended release duration: 0.2

Adjust the suction intensity according to individual endurance.

Recommended application:

- a.Apply essential oil evenly to the treatment area.
- b.Lift slowly from the lower arm towards the elbow, and from the elbow towards the shoulder using the handle.
- c.Slowly move the handle in circular motions from the elbow towards the shoulder.
- d.Slowly lift the upper arm from the inside outward.
- e.Move the handle in circles on the upper arm.
- f.Repeat on the other side.
- g.Wipe arms with a towel.

5.LIPO LASER - Excess adipose tissue elimination and result consolidation

Auxiliary tools: strap

Recommended working time: 20 minutes.

Recommended energy: 3

Recommended Mode: M1

a.Fix the laser pads on the arms, 20 minutes.

b.Remove the laser pads and the strap.

6.Treatment schedule

Following the initial course of 12 sessions, adipose tissue reduction initiates, and the skin gains firmness and tightness. The intermediate treatment spanning 24 sessions contributes to body shaping, clears channels, and enhances overall physical fitness. The advanced treatment, comprising 48 sessions, consolidates and sustains the results, preventing any rebound effects.

Limb Elegance

Handles: CAVI 50K&RF, VACUUM&RF, and Laser Pad

Duration (100 minutes)

Preparation: hot towel, essential oil, gel, strap.

1.Gently wipe the legs with a hot towel.

2.Manual massage - clear channels, enhance circulation

Product required: essential oil

Recommended working time: 10 minutes

Recommended application:

- a.Apply essential oil evenly to the treatment area.
- b.Have the client lie on her stomach and begin with the left leg: Take some oil and warm it up by rubbing it between your hands. Then apply it from the lower leg to the thigh, and down to the heel (soothing act).
- c.Massage the leg from the lower to the upper part using alternating palms, then return to the foot.
- d.Massage the entire leg from lower to upper using alternating first-web strokes with both hands.
- e.Massage the leg in 'S' motions from bottom to top using both hands.
- f.Push 3 meridians from bottom to top with bent fingers.
- g.Soothing act.
- h.Repeat on the other side.
- i.Wipe the leg with a towel.

3.CAVI 50K&RF - clear channels and contour the body

Product required: Gel

Parameters setting:

Recommended working time: 15 minutes

Recommended RF energy: 3

Recommended 50K energy level: 3

Recommended 50K mode: M1

Recommended application:

- a.Apply gel evenly to the treatment area.
- b.Push from the ankle towards the popliteal fossa with the handle.
- c.Move the handle in small circles on the lower leg.
- d.With the RF handle, push from the popliteal fossa to the end of the thigh.
- e.Move the handle in small circles on the thigh.
- f.Lift from the inner side of the thigh outward using your hand and the handle.
- g.Repeat the same on the other leg.
- h.Wipe the legs with a towel.

4.VACUUM&RF - Clarify channels, tighten skin, and sculpt the body

Product required: essential oil

Parameters setting:

Recommended working time: 15 minutes

Recommended RF energy: 3

Recommended suction duration: 0.5

Recommended release duration: 0.2

Adjust the suction intensity according to individual endurance.

Recommended application:

- a. Apply essential oil evenly to the treatment area.
- b. Push from the ankle towards the popliteal fossa with the handle.
- c. Move the handle in small circles on the lower leg.
- d. With the RF handle, push from the popliteal fossa to the end of the thigh.
- e. Move the handle in small circles on the thigh.
- f. Lift from the inner side of the thigh outward using your hand and the handle.
- g. Repeat the same on the other leg.
- h. Wipe the legs with a towel.

5. Manual massage - clear channels, enhance circulation

Product required: essential oil

Recommended working time: 10 minutes

Recommended application:

- a. Let the client lie on his back: Place some oil in your hands, warm it up by rubbing, and then stroke from the foot to the root of the thigh (soothing act).
- b. Alternately push the leg from lower to upper with both palms and then return to the foot.
- c. Massage the entire leg from lower to upper using alternating first-web strokes with both hands.

d. Massage the leg in 'S' motions from bottom to top using both hands.

e. Alternately push the leg from lower to upper using the bent fingers of both hands.

f. Soothing act.

g. Repeat the same on the other leg.

h. Wipe the legs with a towel.

6. CAVI 50K&RF - clear channels and contour the body

Product required: Gel

Parameters setting:

Recommended working time: 15 minutes

Recommended RF energy: 3

Recommended 50K energy level: 3

Recommended 50K mode: M1

Recommended application:

a. Apply gel evenly to the treatment area.

b. Move the handle in small circles on the lower leg.

c. Push from the knee to the root of the thigh.

d. Between the knee and the root of the thigh, move the handle in circles around the leg.

e. Move the handle in small circles on the thigh.

f. Lift from the inner side of the thigh to the middle using both your hand and the handle.

g.Repeat the same on the other leg.

h.Wipe the legs with a towel.

7.VACUUM&RF - Clarify channels, tighten skin, and sculpt the body

Product required: essential oil

Parameters setting:

Recommended working time: 15 minutes

Recommended RF energy: 3

Recommended suction duration: 0.5

Recommended release duration: 0.2

Adjust the suction intensity according to individual endurance.

Recommended application:

a.Apply essential oil evenly to the treatment area.

b.Move the handle in small circles on the lower leg.

c.Push from the knee to the root of the thigh.

d.Between the knee and the root of the thigh, move the handle in circles around the leg.

e.Move the handle in small circles on the thigh.

f.Lift from the inner side of the thigh to the middle using both your hand and the handle.

g.Repeat the same on the other leg.

h.Wipe the legs with a towel.

8.LIPO LASER - Excess adipose tissue elimination and result consolidation

Auxiliary tools: strap

Recommended working time: 20 minutes.

Recommended energy: 3

Recommended Mode: M1

a.Fix the laser pads on the legs, 20 minutes.

b.Remove the laser pads and the strap.

9.Treatment schedule

The primary treatment consists of 12 sessions, focusing on sculpting the legs; the intermediate treatment includes 24 sessions, resulting in shapelier and firmer legs with more noticeable effects; the advanced treatment with 48 sessions not only consolidates the results but also promotes collagen regeneration, making the treatment even more effective.

Dorsal Trimming

Handle: Vacuum&RF

Duration (30 minutes)

Preparation: Hot towel, essential oil.

1.Gently wipe the back with a hot towel

2.Manual massage - clear channels, enhance circulation

Product required: essential oil

Recommended working time: 15 minutes

Massage techniques:

- a.Take some oil and warm it up by rubbing it between your hands. Next, use both hands to stroke the back towards the waist, then return to the neck. Finally, tap the Fengchi and Fengfu points.
- b.Push the Du channel from the neck to the waist side and return to the neck with both hands, then tap the Fengchi and Fengfu points.
- c.Push the Bladder Meridian from the shoulder to the waist side and return to the shoulder with both hands, then tap the Fengchi and Fengfu points (left side first).
- d.Alternately push the whole scapulothoracic joint using the first web of both hands. Then, push it back along the Daban Meridian, ending at the palm, starting with the left side.
- e.Push Daban Meridian from the back of the neck with bent fingers of both hands.
- f.Push Daban meridian using the first web of one hand while stroking it with the other hand, beginning with the left side.
- g.Push from the shoulder to waist side section by section with bent fingers of both hands, then return to the shoulder, and push to the palm.
- h.On one side of the back: Horizontally push Bladder Meridian with thumbs, from the scapula to the waist.
- i.Alternately stroke from the waist side to the shoulder using palms, and end at the palm.
- j.Repeat on the other side.
- k.Wipe the back with a towel.

3.VACUUM&RF - Clarify channels, tighten skin, and sculpt the body

Product required: essential oil

Parameters setting:

Recommended working time: 15 minutes

Recommended RF energy: 3

Recommended suction duration: 0.5

Recommended release duration: 0.2

Adjust the suction intensity according to individual endurance.

Recommended application:

- a.Apply essential oil evenly to the treatment area.
- b.Push the Du channel in 3 sections with the handle.
- c.Push the Du channel with the handle.
- d.Divide the Bladder meridian into 3 sections, then push every section back and forth with the handle.
- e.Push the Bladder meridian with the handle.
- f.Move the handle in circular motions along Bladder meridian to the Baliao point.
- g.Move the handle on the Daban area (the area that connects the neck and shoulders) from the neck.
- h.Lift the scapula seam back and forth.
- i.Move the handle from the left waist side to the right waist side back and forth.
- j.Wipe the back with a towel.

4.Treatment schedule

Upon completion of the initial course (12 sessions), experience relaxation in the shoulders and back, reduced adipose tissue, and relief from the neck hump. The intermediate treatment (24 sessions) focuses on refining the neck hump, clearing meridians, regulating sub-health, and improving sleep. The advanced treatment, spanning 48 sessions, results in a thinner back and more aesthetically pleasing contour, while also enhancing the function of internal organs, promoting digestion, detoxification, and metabolism.

Abdominal Toning

Handles: CAVI 30K&RF, VACUUM&RF, and Laser Pad

Duration (60 minutes)

Preparation: hot towel, essential oil, gel, strap.

1.Gently wipe the abdomen with a hot towel.

2.Manual massage - clear channels, enhance circulation

Product required: essential oil

Recommended working time: 10 minutes

Recommended application:

- a.Apply some oil to your hands and warm it by rubbing them together, then spread it on the abdomen.
- b.Alternately knead the belly with both hands.
- c.Lift from the waist side to the groin with both of your hands alternately.
- d.Repeat on the other side.

- e.Alternately push the ascending colon up, the transverse colon sideways, and the descending colon down using both hands.
- f.Massage around the navel in a clockwise direction with overlapped hands.
- g.Lift the Belt Vessel on both sides of the waist with both hands alternately.
- h.Gently wipe the abdomen with a towel.

3.CAVI 30K&RF - clear channels and contour the body

Product required: Gel

Parameters setting:

Recommended working time: 15 minutes

Recommended RF energy: 3

Recommended 30K energy: 3

Recommended 30K mode: M1

Recommended application:

- a.Apply gel evenly to the treatment area.
- b.Glide the handle in small circular motions around the navel area.
- c.Glide the handle in big circular motions around the navel area.
- d.Move the handle up and down on the waist.
- e.Lift the Belt Vessel on both sides using the handle.
- f.Lift the Belt Vessel on both sides with the handle and your hand.
- g.Gently wipe the abdomen with a towel.

4.VACUUM&RF - Clarify channels, tighten skin, and sculpt the body

Product required: essential oil

Parameters setting:

Recommended working time: 15 minutes

Recommended RF energy: 3

Recommended suction duration: 0.5

Recommended release duration: 0.2

Adjust the suction intensity according to individual endurance.

Recommended application:

- a.Apply gel evenly to the treatment area.
- b.Glide the handle in small circular motions around the navel area.
- c.Glide the handle in big circular motions around the navel area.
- d.Move the handle up and down on the waist.
- e.Lift the Belt Vessel on both sides using the handle.
- f.Lift the Belt Vessel on both sides with the handle and your hand.
- g.Gently wipe the abdomen with a towel.

5.LIPO LASER - Excess adipose tissue elimination and result consolidation

Auxiliary tools: strap

Recommended working time: 20 minutes

Recommended energy: 3

Recommended Mode: M1

- a.Fix the laser pads on the abdomen, 20 minutes.
- b.Remove the laser pads and the strap.

6.Treatment schedule

The initial treatment of 12 sessions enhances the definition of the waistline. The intermediate treatment, spanning 24 sessions, not only sculpts the body but also imparts firmness and plumpness to the skin. The advanced treatment of 48 sessions not only solidifies the achieved results but also provides relief from irregular menstruation, a cold palace, and painful menstruation.

Booty Sculpt

Handle: Vacuum&RF, Laser Pad

Duration (40 minutes)

Preparation: Hot towel, essential oil, strap.

1.Gently wipe the buttocks with a hot towel.

2.Manual massage - clear channels, enhance circulation

Product required: essential oil

Recommended working time: 10 minutes

Recommended application:

- a.Stand beside the client and take some oil and warm it by rubbing between your hands, then apply it from the Baliao Area to the waist, and back to the buttocks (soothing act).
- b.Push Baliao Area with thumbs.
- c.Soothing act.
- d.Left side first, push the entire buttock from the root of the thigh upward with both hands.
- e.Overlap hands and push from the root of the thigh to Belt Vessel.
- f.Massage the buttocks from the side towards the highest point, 3 times (sculpting).
- g.Soothing act.
- h.Repeat on the other side.
- i.Wipe buttocks with a towel.

3.VACUUM&RF - Clarify channels, tighten skin, and sculpt the body

Product required: essential oil

Parameters setting:

Recommended working time: 10 minutes

Recommended RF energy: 3

Recommended suction duration: 0.5

Recommended release duration: 0.2

Adjust the suction intensity according to individual endurance.

Recommended application:

- a.Apply essential oil evenly to the treatment area.
- b.Lift from the root of the thigh towards Belt Vessel with the handle.
- c.Lift the buttocks from the side towards the highest point using the handle.
- d.Move the handle in small circles on the entire buttocks.
- e.Lift from the root of the thigh towards Belt Vessel with the handle.
- f.Lift the buttocks from the side towards the highest point using the handle.
- g.Repeat on the other side.

4.LIPO LASER - Excess adipose tissue elimination and result consolidation

Auxiliary tools: strap

Recommended working time: 20 minutes

Recommended energy: 3

Recommended Mode: M1

- a.Fix the laser pads on the buttocks, 20 minutes.
- b.Remove the laser pads and the strap.

5.Treatment schedule

After completing the initial 12-session regimen, you'll experience enhanced definition in your buttock contours; progressing to the intermediate course with 24 sessions refines the buttocks and promotes skin elasticity. Advancing to the comprehensive 48-session program not only alleviates symptoms related to cold body and cold uterus but also addresses gynecological concerns. It reinstates regular menstruation, solidifies your progress, and elevates your overall feminine allure.

Contact Us



 Aristormofficial  