

LED red light physiotherapy belt

Product introduction

This product is a professional substitute of outdated and bulky medical equipment. It is easier for people with fat to get slim by stimulating the lymphatic system to discharge fat waste from the body. Red light can stimulate the generation of collagen and soften fine wrinkle and plication. It is superior to any other face cream or essence.

After using the product, aerobic exercise of 20minutes is suggested (most importantly, in this step, your body will exclude fat from your lymphatic system and liver). If you have enough time and energy, you can do a few minutes of weight-bearing exercise (such as squatting, push-ups, sit-ups) to maximize the effect in the therapy area.

We strive to solve your problem instead of providing medical advice. The best suggestion on making full use of the product is a physiotherapy for you.

Feel free to contact us: 404126844@qq.com

Usage time:

Slimming: 20 minutes/time, 3-5 times a week; skin rejuvenation: 10 minutes/time, 3-5 times a week; pain-relieving: 20-30 minutes /time, 7 times a week

Chronic fatigue: 10-20 minutes / time, 3-7 times a week; arthritis: 10 minutes / time, 7 times a week; would healing: 10 minutes/time, 3 times a week

Product parameters:

Product size: 160*35cm

Adapter: Input voltage AC 100-240V, use voltage DC 12V 10A

Power control: P1-P5, a total of 5 positions can adjust the power. Start button: long press, switch the machine; short press, frequency switching, 2 levels (10Hz/ 40Hz)


Time control: 10--90 minutes (the best use time is 30 minutes)

Main function: weight loss, relieve soreness, skin rejuvenation

Products Show



布局图样



Button distinction description:

"On/Off key": long press for 3 seconds to switch on/off / short press to switch frequency 10HZ, 40HZ

"Power key": power is divided into 5 grades, P1, P2, P3, P4, P5

"Time key": time can be set, 10 20 30 40 50 60 70 80 90 minutes, 9 time periods can be adjusted