

Model No.: MS-13H1



Positioning: A breast enlargement device

Packing List

- 1 x Main Machine
- 1 x Filter
- 1 x Tube for Cups
- 2 x Pairs of Pink Breast Vacuum Cups
- 2 x Connecting Cord for Cups

Technical Parameters

Max vacuum suction pressure: 70KPa

Microcurrent frequency: 100KHz

Motor rotation speed: 12000 (rpm)


Rated power: 48W

Precautions on Use

1. You must use a plug with a ground pin and ensure a power socket is already grounded before using the equipment.
2. Make sure the equipment's voltage is adaptive. If the local voltage is unstable, we suggest users add a voltage regulator with matching power between the main supply and the equipment.
3. Please uniformly use specified parts provided or indicated by its original manufacturer to guarantee the equipment's curative effect and expected service life.
4. Please do not put the equipment near a water source, in a damp place, or in sites that expose it to sunlight directly.
5. Please do not place the equipment near an intense heat source since this may affect its service life and everyday use.
6. Those with a prosthesis inside their breasts can not use it.
7. Those with implanted metal stents (not including a contraceptive ring placed inside a women's uterus) can not use it.
8. Women in pregnancy or breastfeeding can not use it.
9. You can not use it after drinking alcohol.
10. Those with heavy menstrual blood can not use it.
11. Those who did cesarean in the past six months can not use it.
12. Those who experienced a miscarriage in the past three months can not use it.
13. Those who had natural labor in the past two months can not use it.
14. If a vacuum cup's suction is too strong, please stop using it immediately and reset the suction.
15. After an operation, clean the equipment with normal saline to ensure its cleanliness and hygiene to extend its service life.

General function introduction:


The following initial interface will show up after turning on the equipment.



Return



Start/ Pause



Time Display & Setting (Time range: 0-60 minutes; suggested operating time: 15-20 minutes)



Time Setting (decrease)



Time Setting (increase)




Energy Level Display & Setting for BIO (Energy level range: 0-9; advised energy level: 3-7)



Energy Level Setting (decrease)



Energy Level Setting (increase)



Mode Selection (BIO)

M1: Continuous working. Namely, when clicking Start, the handpiece releases energy all the time.

M2: Discontinuous working. Namely, when clicking Start, the handpiece releases energy on and off.



Start/ Pause for Vibration

ON: Start


OFF: Pause



Vacuum Meter



Rotating knob for setting suction (rotate to the left to decrease suction and to the right to increase)



Suction Time Display & Setting (Time range 0.0-2.0 seconds; suggested operating time: 0.3-1.5 seconds)



Release Time Display & Setting (Time range 0.0- 2.0 seconds; suggested operating time: 0.0-0.8 seconds)



Suction & Release Time Setting (decrease)



Suction & Release Time Setting (increase)



Mode Selection (Vacuum)



Direct suction (select it if you want to do cupping for the long term or to raise transitory petechiae and ecchymosis due to heavy dampness inside your body)



One suction followed by one release (choose it if you want to raise transitory petechiae and ecchymosis or to do short-term cupping)



One suction followed by two releases (select it if you want to massage and to raise transitory petechiae and ecchymosis slightly)



One suction followed by three releases (select it if you want to massage, ease skin, and alleviate soreness)

technical parameter

Maximum pressure: 70KPa
Microcurrent frequency: 100KHz
Vibration motor speed: 12000 (rpm)
Rated power: 48W

Working Principle

Vacuum Cup for Breast Enlargement

It uses low currents to form a magnetic field inside a human body. A beautician doesn't need to use any strength but his hands to enlarge breasts. When the magnetic field passes by the breasts, the breasts will feel slightly numb. It's because the magnetic field deeply dredges and breaks up a blocked part, which activates mammary gland factors and makes them revive, regenerate, and regrow. Meanwhile, it dredges the breasts' channels and collaterals and lymphatic ducts and promotes lymphatic circulation to revive and regenerate breasts' adipose tissue and accumulate and thicken the breasts' fat so that breasts can become larger, fuller, and upstanding.

Principle of breast enlargement with rhythmic physical absorption

1. Based on the principle of physical absorption, the device is developed per human physiology and is the most direct and effective equipment to exercise breasts.
2. It utilizes bionic-physical technology to channel estradiol secreted at the peak of the human body and help breasts become a target organ of estradiol, which activates fat stored in plasma and transfers it to breasts to enlarge their size.
3. The violent micro-movements fully exercise muscular tissue and stimulate the nervous system to secrete estradiol to enlarge the breasts naturally.
4. The device's absorption can exercise the elasticity of muscular tissue to lift the breasts and make them full and upstanding. The suction and release and movements created can also constantly stimulate the pectoralis major and mammary gland.
5. Meanwhile, it can make the breasts accumulate fat instantly.

Principle of breast enlargement with vibration

Micro-movements activate the root of the mammary gland instantly, making under-developed, atrophic, or saggy breasts generate natural movement, causing flat, small, saggy breasts to redevelop, and reshaping upstanding and full breasts.

Indications

1. With slack or saggy breasts.
2. With flat or asymmetrical breasts or whose breasts have outward expansion.
3. Who have nipple retraction and pursue more perfect and healthier breasts.
4. With breast nodules, slight hyperplasia, or who experience distending pain in breasts during menstruation.
5. Who think she has a less developed mammary gland and mastatroph after childbirth.
6. With sunken, flat, and less than sexy hips.

7. Whose hips are flat, shapeless, or with less fat.
8. Who want to have a perfect figure due to psychological or professional needs.
9. Whose hips are in good size but with inadequate plumpness or lack hip peak.
10. With square hips that caused by fat accumulation.
11. With shoulder or back soreness or stiff neck.
12. With insomnia, dreaminess, or a fading memory.
13. Prone to fatigue or drowsiness or whose vital energy and blood circulation are obstructed.
14. With a thick back or who look unsightly in clothes.
15. With Dowager's Hump.
16. With cold hands, cold feet, a cold uterus, or a cold-natured body.
17. With flab on waist and abdomen, potbelly, or flabby skin after delivery.
18. With unsightly waistline due to sitting too long.
19. With striae distensae or stretch marks.
20. With constipation or whose abdomens have obstructed channels and collaterals.

Contraindications

1. Women in menstruation, pregnancy, or breastfeeding.
2. With surgical wounds or in surgical recovery.
3. With heart disease, hypertension, etc.
4. With skin disease, infectious disease, or whose skin is in a sensitive period.
5. Who have a weak body.
6. With metal implants inside the body (such as a stent, pacemaker, etc.) or allergic to metals.
7. With malignant tumors, hemophilia, or severe bleeding.
8. Who are drunk, thirsty, overworked, or with a full or empty stomach.

Healing Crisis:

Numbness--obstructed local circulation

Pricking--inflammation or an old illness

Ache--over one year of obstruction

Sore and swollen--blocked acupoints or muscle synechia

Slight redness--poor subcutaneous circulation

Swelling pain--severe cold-dampness, poor metabolic cycle, or lactic acid metabolizes outward

Erythema--endogenous heat

A skin surface feels painful--dry skin or slack anadesma

Itching--severe cold-dampness

Body parts can use the device: breasts, buttocks, abdomen, and back.

Programs: breast enlargement/ buttocks enhancement/ abdomen sculpting/ back sculpting

Entertaining clients: To know a client's concerns and psychological needs. Listen to his experience, and then recommend a proper program to him.

Case 1

- A. I want to have breast enlargement. Do you have some suitable programs for recommendation?
- B. OK. Let me lead you to a room first. Please change your clothes. I'll offer you the appropriate program after checking your breast conditions.

Preparation:

Towel, essential oil, gel

1. Clean breasts with the towel.
2. Give proper advice per a client's breast conditions.
e.g., Your breast bottom looks nice, but it has too much free fatty acid. I advise you to try a "breast enlargement" program. It takes 40 minutes. I'll massage you to make you relaxed and promote blood circulation. Then I'll place vacuum cups for breast enlargement on your breasts. The cups unite vacuum and BIO, which can better dredge channels and collaterals and gather the breasts' fat. Just enjoy yourself and tell me immediately if you are not comfortable with the suction level.

Procedures:

Now I am going to start, and the main functions of this treatment are:

- Relieve breast with loosening, sagging, or outward expansion.
 - Alleviate distending pain in breasts before menstruation.
 - Dredge and adjust breast.
 - Prevent galactoceles.
 - Prevent breast cancer.
 - Breast Sculpting
1. **Measure the size and take photos.**
 2. **Massage** (dredge channels and collaterals and make the client feel eased and relaxed)

Techniques:

- a. Stand next to the head of a bed, apply essential oil from Danzhong (RN17) to an armpit with hands, and lift Cooper's ligaments while moving the hands upward. Repeat three times.
- b. Press Danzhong (RN17), Ruzhong (ST17), Dabao (SP21), Yinchuang (ST16), Zhongfu (LU1), and Yunmen (LU2) with thumbs. Repeat three times.
- c. Caress the treated part. Repeat three times.
- d. Draw Arabic numeral eight between Danzhong (RN17) and breasts with palms overlapped. Repeat three times.
- e. Lift from accessory breasts to Cooper's ligaments with hands; do it in left-right order. Repeat ten times.
- f. Caress the treated part. Repeat three times.
- g. Sit next to the client; caress and lift her breasts with hands; do it in left-right order. Repeat 3 to 5 times.

- h. Push lactiferous ducts around the breasts with hands' purlicue. Repeat 3 to 5 times.
- i. Dredge the breasts' nodules with thumbs moving circlewise and alternately. Repeat three times.
- j. Caress and lift the whole breast with hands. Repeat 3 to 5 times.
- k. Do the same on the other side.

3. Vacuum Cup for Breast Enlargement (dredge and correct breasts)

Advised operating time: 15-20 minutes; suggested energy level for BIO: 3-7, and recommended mode for BIO: M1. Vibration Mode: ON. Advised suction time: 0.3-1.5 seconds; suggested release time: 0.0-0.8 seconds; recommended mode: one suction followed by one release; suction can be adjusted per the client's tolerance.

- a. Vacuum cup for breast enlargement: select proper cups for the client, and fasten them onto her breasts for around 20 minutes.
- b. Remove the cups.
- c. Clean the treated part with a towel. Treatment ends.

Remark: The vacuum cups can be customized per the breast size of the client.

The course of treatment for guidance: to ask the client about her feelings after treatment, whether she wants to take it for the long term, and her expectations; then customize a private course of treatment per her needs. When you accomplish a treatment cycle, breasts will be lifted. After a course, or 12 times of treatment, the breasts become firm and enlarged; when you finish 24 times of treatment, accessory breasts, outward expansion, and sagging start to tighten and concentrate slowly. If you complete 48 times of treatment, the breasts turn large and upstanding. If you stick to the treatment, the effect will get more noticeable.

Matters Needing Attention After Treatment

- 1. Take a shower after 4 to 6 hours.
- 2. Avoid exposing yourself to a windy environment and catching a cold.
- 3. Avoid eating spicy or greasy food.
- 4. Drink more warm water to expedite metabolism.
- 5. Avoid staying up late, smoking, or drinking.

Case 2

- A. Do you have programs for buttocks lifting and tightening? I want to have a try.
- B. OK. Let me lead you to a room first. Please change your clothes. I'll offer you the appropriate program after checking your buttocks conditions.

Preparation:

Towel, essential oil, gel

- 1. Clean buttocks with the towel.
 - 2. Offer sound advice per a client's buttocks conditions.
- e.g., Your buttocks' fat show outward expansion. I advise you to try "buttocks enhancement." It takes 40 minutes. I'll massage you to dredge channels and

collaterals and promote blood circulation. Then I'll place vacuum cups for breast enlargement on your breasts. The cups unite vacuum and BIO, which can better dredge channels and collaterals and gather the buttocks' fat. You just enjoy yourself and tell me immediately if you are not comfortable with the suction.

Procedures:

Now I am going to start, and the main functions of this treatment are:

- Relax muscle and boost local blood circulation.
- Prevent amyotrophy and relieve muscle spasms.
- Firm skin and build the body.
- Relieve sagging caused by slack muscles.
- Stimulate muscle contraction, make muscles plump and muscular, and build sexy buttocks.

1. Measure the size and take photos.

2. Massage (dredge channels and collaterals and make the client feel eased and relaxed)

Techniques:

- a. Stand sideways with hands applying essential oil, slide from the sacral region (BL31-BL34) to the waist, and lift upwards along buttocks. Repeat three times.
- b. Push the sacral region (BL31-BL34) with thumbs. Repeat three times.
- c. Caress the treated part three times, and then press Shenshu (BL23), sacral region (BL31-BL34), Changqiang (DU1), Huantiao (GB30), and Chengfu (BL36). Repeat three times.
- d. Caress the treated part. Repeat three times.
- e. Do it in left-right order; push from thigh root to Daimai (GB26) with hands along Bladder Meridian, Kidney Meridian, Liver Meridian, and Gallbladder Meridian. Repeat three times, respectively.
- f. Push upwards from the thigh root to Daimai (GB26) with hands along Bladder Meridian, Kidney Meridian, Liver Meridian, and Gallbladder Meridian. Repeat three times, respectively.
- g. Push from the thigh root to Daimai (GB26) with palms overlapped (lifting). Repeat five times.
- h. Push one line after another and back and forth from the buttocks' two sides to their highest point with hands (shaping). Repeat three times.
- i. Repeat technique No.7.
- j. Caress the treated parts.
- k. Do the same on the other side.
- l. Rub the sacral region (BL31-BL34) with your hands till it turns hot.

3. Vacuum Cup for Breast Enlargement (dredge and lift buttocks)

Advised operating time: 15-20 minutes; suggested energy level for BIO: 3-7, and recommended mode for BIO: M1. Vibration Mode: ON. Advised suction time: 0.3-1.5 seconds; suggested release time: 0.0-0.8 seconds; recommended mode: one suction followed by one release; suction can be adjusted per the client's tolerance.

- a. Do moving cupping with vacuum cup for breast enlargement: lift from a thigh

- root to the highest point of the buttocks. Repeat 3 to 5 times.
- b. Lift from the buttocks' two sides to the middle. Repeat 3 to 5 times.
 - c. Fixed cupping: Select proper cups, fasten them onto buttocks, set the energy to a reasonable level, and leave them there for around 5 to 15 minutes.
 - d. Remove the cups and clean the treated part, two minutes.

Remark: The cup can be customized per the breast size of the client.

4. Remeasure the breasts, take photos, and compare them.

The course of treatment for guidance: to ask the client about her feelings after treatment, whether she wants to take it for the long term, and her expectations; then customize a private course of treatment per her needs. When you accomplish a treatment cycle, buttocks will be lifted. After a course, or 12 times of treatment, excess flab on the buttocks will vanish; when you finish 24 times of treatment, fat around the buttocks start to concentrate. If you complete 48 times of treatment, your buttock curve will form, their elasticity will increase, and your gynecological diseases will be relieved.

Matters Needing Attention After Treatment

1. Take a shower after 4 to 6 hours.
2. Avoid exposing yourself to a windy environment and catching a cold.
3. Avoid eating spicy or greasy food.
4. Drink more warm water to expedite metabolism.
5. Avoid staying up late, smoking, or drinking.
6. It's advised to wear hip-lifting pants to prevent loosening and sagging.
7. Do some simple exercise to lift your buttocks.

Case 3

- A. A: I want to have abdomen sculpting. Do you have some programs for recommendation?
- B. B: OK. Let me lead you to a room first. Please change your clothes. I'll offer you the appropriate program after checking your abdomen conditions.

Preparation:

Towel, essential oil, gel

1. Clean abdomen with the towel.
2. Give some good advice per the client's abdomen condition.
e.g., You have a potbelly, and your skin is slack. I advise you to try "abdomen sculpting." The program takes 40 minutes. I'll massage you to dredge channels and collaterals and promote blood circulation. Then I'll do cupping with vacuum cups for breast enlargement. It works well on curing a common cold caused by wind-cold and a cold uterus and can regulate the body's internal organs, balance Yin and Yang, and enhance immunity. Just enjoy yourself and tell me immediately if you are not comfortable with the suction.

Procedures:

Now I am going to start, and the main functions of this treatment are:

- Relieve a woman's cold hands, cold feet, cold uterus, or cold-natured body.
- Relieve waist and abdomen's flabby and soft skin.
- Remove flab from waist and abdomen, and relieve swim-ring-like belly and potbelly.
- Tighten the skin, reduce stretch marks and striae distensae, and increase skin elasticity.
- Accelerate metabolism, relieve constipation, and enhance the movement of the intestinal tract.

1. Measure the size and take photos.

2. **Massage** (dredge channels and collaterals and make the client feel eased and relaxed)

Techniques:

- a. Rub essential oil into the abdomen with hands moving circlewise. Repeat 3-5 times.
- b. Rub the belly back and forth with your hands. Repeat 3 to 5 times.
- c. Rub the surrounding area of the belly button with combined hands drawing small circles. Repeat 16 times.
- d. Rub the surrounding area of the belly button with combined hands drawing big circles. Repeat 16 times.
- e. Stroke the two sides below the waist towards ribs. Repeat 3-5 times.
- f. Stroke the two sides below the waist towards groin. Repeat 3-5 times.
- g. Push from the side of the waist to the belly button with hands doing it alternately. Repeat 3 to 5 times.
- h. Push from the side of the waist to the groin with hands doing it alternately. Repeat 3 to 5 times.
- i. Rub Daimai (GB26) of the two sides as the way of rubbing the spine with hands doing it alternately; do it in left-right order. Repeat 10-15 times.

3. **Vacuum Cup for Breast Enlargement** (dredge and tighten)

Advised operating time: 15-20 minutes; suggested energy level for BIO: 3-7, and recommended mode for BIO: M1. Vibration Mode: ON. Advised suction time: 0.3-1.5 seconds; suggested release time: 0.0-0.8 seconds; recommended mode: direct suction; suction can be adjusted per the client's tolerance.

- a. Fixed cupping: Select proper cups, fasten them onto the abdomen, and leave them there for around 5 to 15 minutes.
- b. Remove the cups and clean the treated parts, two minutes.

4. **Remeasure the size, take photos, and compare them.**

The course of treatment for guidance: to ask the client about her feelings after treatment, whether she wants to take it for the long term, and her expectations; then customize a private course of treatment per her needs. When you accomplish a treatment cycle, the abdomen is heated, which expedites metabolism. After a course, or 12 times of treatment, fat metabolizes slowly, and constipation is relieved. When

you complete 24 times of treatments, the effect becomes more apparent, and the skin turns tightened and is lifted. If you finish 48 times of treatment, excess fat disappears gradually, the waist and abdomen's curves sculpt, and a slim waist forms.

Matters Needing Attention After Treatment

1. Keep the abdomen warm and avoid exposing to a windy environment and catching a cold.
2. Avoid eating and drinking too much, staying up late, drinking alcohol, or eating raw, cold, spicy, or greasy food. Instead, drink more hot water.
3. Take a shower after 4 to 6 hours.
4. Avoid enjoying sauna or hot springs or doing strenuous exercise within seven days.
5. Rub clockwise your abdomen with your hands before sleep at night, achieving better weight-loss and metabolism effect.

Case 4

- A. I want to have back sculpting. Do you have some suitable programs for recommendation?
- B. OK. Let me lead you to a room first. Please change your clothes. I'll offer you the appropriate program after checking your back conditions.

Preparation:

Towel, essential oil, gel

1. Clean back with the towel.
2. Offer sound advice per a client's back conditions.
e.g., You have a thick back and heavy dampness. I advise you to try "back sculpting." The program takes 40 minutes. I'll massage you to dredge channels and collaterals and promote blood circulation. Then I'll do cupping with vacuum cups for breast enlargement. It works well on curing a common cold due to wind-cold and shoulder and neck strain and can regulate the body's internal organs, balance Yin and Yang, and enhance immunity. You just enjoy yourself and tell me immediately if you are not comfortable with the suction.

Procedures:

Now I am going to start, and the main functions of this treatment are:

- Alleviate shoulder and back soreness and relieve Dowager's Hump.
- Dredge channels and collaterals and relieve clogged channels and collaterals.
- Accelerate blood circulation and metabolism.
- Improve head blood supply and sleep.
- Regulate viscera functions and strengthen the physique.
- Tighten skin and prevent slack and soft skin.
- Remove excess flab from the back and sculpt the back.
- Improve metabolism.

- Promote the immune system's functions.

1. Measure the size and take photos.

2. **Massage** (dredge channels and collaterals and make the client feel eased and relaxed)

Techniques:

- Apply oil to the whole back with hands and press Fengchi (GB20) and Fengfu (DU16). Repeat three times.
- Push Du Meridian with hands doing it alternately. Repeat three times.
- Push Bladder Meridian (BL) with hands doing it alternately; do it in left-right order. Repeat three times.
- Push from Dazhui (DU14) to Du Meridian and Bladder Meridian with the base of palms doing it alternately, and lift from the side of waist till the palms slide out from armpit, and do it in left-right order. Repeat three times.
- Push from the back of neck and shoulder inch by inch with kneeling fingers. Repeat 5-8 times.
- Start with the back of the neck, pushing anadesma connecting the neck and shoulder with kneeling fingers. Repeat 5-8 times.
- Push from the back of neck to waist inch by inch with kneeling fingers. Repeat 3-5 times.
- Stroke anadesma connecting the neck and shoulder with a hand's purlicue and caress it with the other hand. Repeat 5-8 times.
- Do the same on the other side.
- One side: push the whole medial border of the scapula with hands' purlicue doing it alternately, repeat 5-8 times; push from the medial border of the scapula to shoulder with hands doing it by turns, repeat 5-8 times; slide out along anadesma connecting the neck and shoulder.
- Caress from the side of the waist to shoulder with palms doing it alternately; press Jianjing (GB21); slide out from arms to palm. Repeat three times.
- Do the same on the other side.

3. Vacuum Cup for Breast Enlargement (dredge and tighten)

Advised operating time: 15-20 minutes; suggested energy level for BIO: 3-7, and recommended mode for BIO: M1. Vibration Mode: ON. Advised suction time: 0.3-1.5 seconds; suggested release time: 0.0-0.8 seconds; recommended mode: one suction followed by one release; suction level can be adjusted per the client's tolerance.

- Push Du Meridian directly. Repeat 3-5 times.
- Push Bladder Meridian directly. Repeat 3-5 times.
- Move from the left side of the waist to the right side, and do it back and forth. Repeat 3-5 times.
- Readjustment: It's direct suction when the release time is set to zero.
- Place vacuum cups on Dazhui(Du-14), Tianzong(SI-11) of the two sides, sacral region (BL31-BL34). Three minutes for each acupoint.
- Remove the cups and clean the back with a towel.

4. Remeasure the size, take photos, and compare them.

The course of treatment for guidance: to ask the client about her feelings after

treatment, whether she wants to take it for the long term, and her expectations; then customize a private course of treatment per her needs. When you accomplish a treatment cycle, the back gets relaxed. After a course, or 12 times of treatment, the back and shoulder become eased, back fat lessens, and Dowager's Hump relieves. When you complete 24 times of treatment, Dowager's Hump improves, channels and collaterals dredges, and the sub-health conditions and sleep take a turn for the better. If you finish 48 times of treatments, your back gets thinner, outlines get nice, viscera function enhances, and digestion, detox, and metabolism boost.

Remark: the cups are too big and only suit clients with thick back. If your back is thin, you can customize the cups.

Matters Needing Attention After Treatment

1. Keep warm, avoid exposing to a windy environment and catching a cold, and drink more hot water.
2. Take a shower after 4 to 6 hours.
3. Avoid staying up late, drinking alcohol, and eating and drinking too much.
4. Avoid eating raw, cold, and spicy food; have enough sleep.
5. Avoid wearing shoulder-baring and backless clothes.

Troubleshooting

1. **I can't start the equipment, and its indicator lamp in the back doesn't light up.**
 - A. Make sure you plug the power cord into a good power socket.
 - B. Check whether the fuse tube in the back is loose or burnt out.
2. **The equipment doesn't have suction, or the suction becomes weak.**
 - A. Check whether the handpiece and the plug connecting to the machine body are closely linked.
 - B. Please check whether the oil filter is tightened and whether the rubber ring is worn out. Air leakage in this place may cause insufficient air pressure.
 - C. If the solution mentioned above doesn't work, please contact the distributor for assistance.
3. **I can start the equipment, but an error message shows on the monitor.**
 - A. Take out the power plug from the back, wait about one minute, reinsert the plug, and restart the equipment.
 - B. If the solution mentioned above doesn't work, please contact the distributor for assistance.

FAQs

1. Q: Does breast enlargement with vacuum cups have side effects?

A: No. It adopts a purely physical method for breast enhancement. Massage and vacuum suction and release are all the commonly used principles of physiotherapeutic

instruments. This equipment, however, applies high-tech and unites various principles, which has an instant effect. It's a breast enhancement device that doesn't have side effects, and its effect will not rebound.

2. Q: How to detox with a vacuum?

A: Vacuum can congest capillaries and stimulate cells to increase their vitality. During operation, the air suction and release and the alternate increase and decrease of the vacuum make local pores open and close continually, promote skin breathing, increase skin oxygen uptake, and speed up the elimination of wastes. Therefore, it can invigorate Qi, activate blood circulation, remove cold and dampness, dredge channels and collaterals, and expel toxins and heat.