User Manual

BS-33D1

Content

1. Safety Precautions	1
2. Contraindications	1
3. Product Overview	2
4. Packing List	10
5. Technical Specification	11
6. Key Advantages	11
7. Operation Guide	11
8. Treatment Schedule	19
9. Use Precautions	20
10. Maintenance & Cleaning	21
11. FAQ	21
12. Ask & Answer	22
13. Warranty & Support	24

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1. Safety Precautions

1. Please read this manual carefully before use and familiarize yourself with the operation of this instrument and related technology.
2. This physiotherapy instrument is a high-precision electronic device that requires proper storage. Please avoid storing it in humid conditions, direct sunlight, or near flammables and explosives to prevent damage.
3. Always use the designated accessory provided or recommended by the manufacturer to ensure performance throughout its service life.
4. Do not place the instrument in high-temperature conditions, as this may negatively impact its service life and performance.
5. Do not apply the instrument on the thyroid, parathyroid, testicles, pregnant women's abdomen, pacemakers, etc.
6. Patients should refrain from using this instrument unless they have obtained permission from their doctor.
7. Ensure that the instrument is applied on moist skin and avoid the operation on dry skin.
8. Power off the instrument when not in use to ensure its safety.
9. If you are currently taking anti-adiposity drugs, please discontinue their use for 1 - 2 months before using this instrument. Failure to do so may result in an extended course of treatment.
10. To maintain hygiene and extend the instrument's service life, sterilize it with saline.

2. Contraindications

Before using this instrument, it is essential to consult your doctor or seek advice from professionals. Details below:

1. Pregnant or lactating woman.
2. Individuals with hypertension, heart disease, or a pacemaker installed.
3. Individuals with unhealed wounds or those in the recovery period from surgery.
4. Individuals suffering from epilepsy, severe diabetes, and hyperthyroidism.
5. Individuals suffering from malignant tumors, hemophilia, and hemorrhage.
6. Individuals with skin disease or infectious disease.
7. Individuals who have had a cesarean section within the last six months, a miscarriage within the last three months, or a natural birth within the last two months should refrain from using this instrument.
8. For your safety, refrain from using if you've had plastic surgery or have implants/metal objects near the treatment area within 3 months.

3. Product Overview

1. Product Introduction

Recognized as a premier physiotherapy instrument, it offers a myriad of advantages, such as alleviating inflammation, providing pain relief, activating the skin, fostering wound healing, and improving circulation. This versatile tool effectively replaces various body and facial solutions typically found in salons. It stands out for its safety and efficacy, eliminating the need for injections, medications, or surgery. The entire process is devoid of pain, invasiveness, or undesirable side effects, ensuring a gradual yet noticeable enhancement through external application. Adopt genuine salon-grade physiotherapy equipment with

native 448KHz frequency, equipped with 2 RET probes and 3 CET probes of different interchangeable sizes, together with fascia blades designed for all body areas.

2. Working principle

CET handle

Capacitive Electrical Transfer Mode: CET mode acts on soft tissues with high water content, such as muscles, fascia, internal organs, lymphatic and vascular systems, with a focus on local and superficial effects.


RET handle

Resistive Electrical Transfer Mode: Designed for hard tissues with lower water content, such as ligaments, tendons, joints, and bones, targeting the deep areas of the body and involving multiple joints

Fascia knife handle


The fascia knife softens the fascia through constant-temperature heating and loosens adhesions with its physical blade head, which can deeply relax muscles, relieve soreness, and promote blood circulation.

3. Comprehensive treatment


RET:

Target tissues: Deep muscles, tendons, ligaments, joint capsules, periosteum
 Target areas: Lower back, buttocks, hips, knees, shoulders, Achilles tendons, stubborn pain points
 Depth: Penetration depth 8-10 cm

CET :


Target Tissues: Subcutaneous fat, superficial fascia, skin, lymph nodes
 Target Areas: Face, neck, edema, superficial muscles, body shaping
 Depth: Penetration depth 3-4 cm

Fascial knife:

Main tissues targeted: Deep fascia, myofascial adhesions, deep muscles, tendons, nodules
 Main areas targeted: Lower back, buttocks, legs, strained fascia
 Depth: Penetration depth 8-10 cm

4. Replace the treatment head


CET and RET: Select the desired treatment head, insert it into the slot of the hand grip and push it firmly into place.

Note: clean the steel tube with a brush or soft cloth before the replacement.

5. Fascia Knife Installation

- 1) Remove the aviation connector from the RET handle on the main unit.
- 2) The fascia knife is connected to a clamp to secure the fascia knife head.
- 3) Connect the fascia knife to the aviation connector on the back of the chassis at the RET position.


6. Interface Introduction


Welcome interface

Upon first use, a language selection interface will appear. You can choose the language according to your needs. Click  to confirm the selected language.


Main Menu

Tap  enter the language selection interface

Tap  and enter the following interface

Main Menu

Reset button

:Suitable for small areas of the body to operate on

Time display and set buttons (0 - 60 minutes)

:Suitable for large areas of the body to be operated on

Energy display and set buttons (0 - 10 levels)

:Click to enter the fascia knife interface

Start/Pause

Tap

and enter the following interface

Main Menu

Reset button

:Eye

Time display and set buttons (0 - 60 minutes)

:Face

Energy display and set buttons (0 - 10 levels)

:Neck

Start/Pause

4. Packing List

No.	Name	Quantity
1	RET handle	1
2	Treatment head of RET handle	2
3	CET handle	1
4	Treatment head of CET handle	3
5	Negative Plate Connection Cable	1
6	Power cord	1
7	Holder	1
8	Pedal Switch	1
9	Fascia Knife	1
10	Fascia knife connecting wire clamp	1

5. Technical Specification

Input Voltage: 110V 60Hz or 220V 50Hz

Rated power: 300W

Frequency: 448Hz

Net weight: 6.3KG

Gross weight: 8.4KG

Product dimensions: 43*32*26CM

6. Key Advantages

1. Achieve superior results in physical therapy by promoting circulation, fostering healing, and encouraging tissue regeneration.
2. Experience inflammation alleviation and pain relief with a safe and effective procedure that avoids discomfort or invasions.
3. Benefit from a unique assembly method for the treatment head, providing a convenient and swift usage experience.

7. Operation Guide

Facial, Eye, and Neck treatment

1. Prepare the instrument, assemble the components, and connect the power.

2. Apply the negative electrode patch to the customer's back.
3. Turn on the power switch on the front of the instrument.
4. Select the CET function in the interface, and then select the handle (choose the appropriate working head).
5. Select the energy level (it is recommended to start from level 3, but adjust according to the level that suits each person).
6. Time suggestion: 10-20 minutes
7. Put on PVC insulating gloves for both hands.
8. Press the "Start" button, place the handle on the control area first, and then press the foot pedal to operate.

Face

Parameters setting

Recommended handle: CET (M)

Recommended working time: 10-20 minutes.

Recommended energy: 3

Product required: massage cream

Eye

Parameters setting

Recommended handle: CET (S)

Recommended working time: 10-20 minutes.

Recommended energy: 3

Product required: massage cream

Neck


Parameters setting

Recommended handle: CET (L)

Recommended working time: 10-20 minutes.

Recommended energy: 3

Product required: massage cream

Body treatment:

1. Prepare the instrument, assemble the components, and connect the power.
2. Interface Functionality Suggestion: RET
3. Apply the negative electrode patch to the customer's abdomen or back.
4. Turn on the power switch on the front of the instrument.
5. Select the energy level (it is recommended to start from level 3, but adjust according to the level that suits each person).
6. Put on PVC insulating gloves for both hands.
7. Time suggestion: 20-30 minutes
8. Press the "Start" button, place the handle on the control area first, and then press the foot pedal to operate.
9. After the RET time is over, the negative electrode plate does not need to be removed. Then select the fascia knife interface.
10. Select the energy level (it is recommended to start from level 3, but adjust according to the level that suits each person).
11. Time suggestion: 10-20 minutes
12. Press the "Start" button, place the handle on the control area first, and then press the foot pedal to operate.

Arm


Parameters setting

1) Recommended handle: RET (M)

Recommended working time: 20 minutes

Recommended energy: 3

Product required: massage cream


2) Recommended handle: Fascia knife

Recommended working time: 20 minutes

Recommended energy: 3

Product required: massage cream

Belly


Parameters setting

1) Recommended handle: RET (L)

Recommended working time: 20 minutes

Recommended energy: 3

Product required: massage cream


2) Recommended handle: Fascia knife

Recommended working time: 20 minutes

Recommended energy: 3

Product required: massage cream

Back


Parameters setting

1) Recommended handle: RET (L)

Recommended working time: 20 minutes

Recommended energy: 3

Product required: massage cream


2) Recommended handle: Fascia knife

Recommended working time: 20 minutes

Recommended energy: 3

Product required: massage cream

Leg


Parameters setting

1) Recommended handle: RET (M)

Recommended working time: 20 minutes

Recommended energy: 3

Product required: massage cream


2) Recommended handle: Fascia knife

Recommended working time: 20 minutes

Recommended energy: 3

Product required: massage cream

8. Treatment Schedule

Treatment area	Treatment duration	Treatment interval	Total sessions
Face	10-20mins	7 days	10 sessions
Eye	10-20mins	7 days	10 sessions
Neck	10-20mins	7 days	10 sessions
Abdomen	20-30mins	2 sessions per week	12 sessions
Arm	20-30mins	2 sessions per week	12 sessions
Back	20-30mins	2 sessions per week	12 sessions
Leg	20-30mins	2 sessions per week	12 sessions

9. Use Precautions

1. Strictly follow the instructions in the manual to use the instrument and train operator.
2. Before using the instrument, use a plug with a grounding pin and ensure the power socket is grounded.
3. Make sure the voltage of the instrument is adapted. If the local power supply voltage is unstable, we recommend adding a power-matched regulated power supply between the mains and the instrument.
4. Remove all the metal objects inside the body to avoid causing accidents or affecting treatment effects.
5. Do not use the instrument on an empty abdomen, and take the treatment 1 hour later after the meal.
6. During operation, the instrument should be in full contact with the skin to avoid uneven heating that could cause burns.
7. Please note that energy output only occurs when you step on the pedal.
8. When using CET/RET, please press the handle firmly against your skin first, and then step on the foot pedal while the handle is moving to avoid burning the customer.
9. The negative electrode plate should be in close contact with the customer's skin. Do not place it through clothing, otherwise there will be no energy output.
10. Do not leave the handle in one place while it is in operation to avoid burning customers.
11. Please adjust the energy level from low to high. During operation, please ask the customer how they feel. If it gets too hot, please lower the energy level immediately.
12. It is recommended that operators wear gloves when operating the instrument to avoid contact with radio frequency current and prevent injury to the operator.
13. Do not drop the handle head on the ground to prevent damage.
14. Do not use the CET handle if the head is damaged to prevent burns.

10. Maintenance & Cleaning

1. Disconnect the power cord when not in use and store the instrument in a dry, cool, and well-ventilated place.
2. Clean and sanitize the applicator with a damp cloth or saline, and then dry it thoroughly with a dry towel to ensure it remains dry.
3. After cleaning the instrument, cover it with a dustcloth or plastic wrap. Before use, ensure to sanitize the handle again for hygienic reasons.
4. Disassembling the instrument and attempting to clean the inside is strictly prohibited. This task should only be performed by professionals.

11. FAQ

1. When can I expect to observe the results?

A: Normally, noticeable results can be observed on the day of the treatment and within a week afterward. Grounded in endogenesis heat therapy, this device employs capacitance and resistance to transmit energy, operating on the internal tissues of the human body. It triggers the natural reparative process, delivering an anti-inflammatory effect. The visibility of results tends to increase with the number of sessions you undergo.

2. Is it safe for the skin?

A: Yes, it is a non-surgical method that stimulates collagen regeneration and metabolism in the deeper layers of the skin without causing any harm. Temporary heat and redness in the treated area are

normal, indicating increased circulation, and should subside within a few hours. There is no need to be concerned about this reaction.

After using the fascia knife handle, some redness, swelling, petechiae, and bruising are normal phenomena and will usually recover within a week.

3. Any side effects?

A: No, there are no side effects. This instrument aims to promote tissue regeneration and circulation while delivering benefits like inflammation alleviation and pain relief.

12. Ask & Answer

Ask	Cause	Answer
1. Fail to turn on the instrument.	→ The power cord is disconnected from the wall socket.	→ Check if the power cord is well connected to the wall socket and ensure that the socket is functional.
	→ The fuse is loose or burned.	→ Replace the fuse on the back of the instrument and restart it.
2. The instrument is not working and can't be restarted.	→ The fuse may have burned out due to unstable voltage or prolonged use.	→ Please replace the fuse and restart the instrument 5 minutes later.

Ask		Cause		Answer
3. The screen shows an error message.	→	Data error.	→	1. Turn off the instrument, disconnect the power cord, and wait for 1 minute to restart the instrument. 2. If the above method doesn't resolve the problem, please contact your dealer.
4. The instrument is not operating.	→	The handle is not properly connected.	→	Please confirm if the handle plug is well connected to the main unit.
	→	Oil in the treatment area could result in insufficient contact between the treatment head and the body.	→	Apply some gels to the treatment area after cleaning.
	→	The treatment head is not in proper contact with the skin.	→	Ensure complete contact between the treatment head and the skin.

13. Warranty & Support

Warranty Policy

1.The main machine is guaranteed for one year and the accessories (like handles and screen) are guaranteed for three months from the date of shipment (by the shipping record).

Within warranty,the seller will deliver free new machine or new accessories if necessary.

Beyond warranty,the seller offers free instruction and troubleshooting service but the repair is with charge (the buyer need to take the cost and the shipping cost for the replacement).

No free service provided for the cases:

- a) Breakdowns caused by personal disassembling and reassembling.
- b) Breakdowns caused by accidental falling or hitting.
- c) Breakdowns caused by a lack of reasonable maintaining.
- d) Breakdowns caused by wrong operation (not following the right guidance of the instruction manual).

2.After purchase, if the machine's parameters or internal structure are modified, we will not provide additional notification services.

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