
NAISIGOO™

IPL Cooling Hair Removal Device **EN**

User Manual

Thank you for choosing this product. In order to facilitate the use, please read this manual carefully and keep it properly.

The pictures in this manual are for illustration only.

CONTENTS

1. Matters Need Attention	3
1.1 Contraindications.....	3
1.2 Warnings and precautions.....	4
1.3 Adverse reaction.....	9
2. Product introduction	9
2.1 Indications for use	9
2.2 About IPL Cooling Hair Removal Device	9
2.3 Packing list.....	10
2.4 Product components.....	11
3. How it works.....	12
3.1 How Does IPL Work.....	12
3.2 Treatment Schedule	12
4 How To Use	12
4.1 Application range	12
4.2 Patch test.....	14
4.3 Operating Steps.....	14
4.4 Skincare after treatment.....	19
5.Troubleshooting	19
6. Technical parameters and specifications.....	22
7. Maintenance	23
8. Disposal.....	23
9. EMC statement.....	24
10. Symbol meaning.....	25
11. Warranty.....	25
12. Contact information	27
Appendix A. EMC information	27

1 Matters Need Attention

Please first read these matters carefully and then follow the user manual to operate. These matters ensure the operation of the product is correct and safe, and prevent users or other people from personal injuries or property damage, please make sure to follow them.

1.1 Contraindications

- Not for individuals under the age of 18. The effects of using the device in this age group are unknown;
- Not suitable for pregnant women, lactating women, women under menstrual periods, individuals with body scars, diabetes, fragile skin and light-sensitive skin;
- Individuals with severe heart disease, moderate to severe hypertension, malignant tumors, folliculitis, infectious skin disease, burned skin and anticoagulant disorders are prohibited;
- Individuals with abnormal collagen metabolism, including history of scars and poor wound healing are prohibited;
- Individuals with vascular diseases are prohibited, such as varicose veins and vasodilation in the irradiation site;
- Individuals with skin sensitive to light and prone to rashes or allergies are prohibited;
- Individuals with irradiated area infected, herpes simplex, skin lesions and hematoma are prohibited;
- Individuals with suppressive immune response diseases (including HIV infection, lupus erythematosus, porphyria are prohibited;
- Individuals with bleeding disorder are prohibited, and those who are using isotretinoin, fruit acid (AHA), salicylic acid (BHA) or azelaic acid are prohibited;

- Individuals taking aspirin drugs, immunosuppressive drugs and other drugs are prohibited. Individuals recently using light-sensitive drugs are prohibited (such as tetracycline and sulfa).

1.2 Warnings and precautions

- Please read the “Instruction Manual” carefully and keep it for future reference;
- If you have any questions during use, please contact the online customer service in the required language;
- The patient is an intended operator;
- The device is intended for a single user;
- This equipment is not intended for use by persons with reduced physical, sensory or mental capabilities, unless they have been given supervision or instruction concerning use of the equipment to avoid the risk of fire or burns;
- Check if the voltage indicated on the equipment corresponds to the local mains voltage before you connect the equipment to avoid risk of electrocution or permanent damage to the equipment;
- If the socket used to power the equipment has poor connections, the plug of the equipment becomes hot. Make sure you plug the equipment into a properly installed socket to avoid the risk of fire and burns;
- Be sure to assemble and use the product in accordance with the instructions contained in the Instruction Manual;
- Do not use accessories, detachable parts and materials not described in the user manual;
- The light exit prolonged contact with the human body may cause irritation, please clean it before and after every use to reduce such potential risk.
- Do not apply any gels or care products (e.g., aloe vera gel) before use;
- Do not use in the bath or shower;

- Do not use this product near water-filled bathtubs, shower basins, washbasins or other water-containing utensils. Please keep this product and power adapter dry;
- If the product is damaged, please do not touch any internal parts to avoid electric shock;
- Do not insert any objects into the product, and do not place the product at will;
- Do not use this product around your eyes to avoid eye injury. Do not use this device to remove eyelashes or eyebrows. Wearing protective glasses is a normal operation that protects the user. If you operate around your eyes, protective glasses may not completely block the energy, and the light will be absorbed by the iris melanin in your eyes at close range;
- We recommend that you wear protective glasses every time you use this product. Although light will not directly affect your eyes, high-frequency flashing light may affect your vision;
- When you need back hair removal, you may need an assistant or partner. This is our recommendation, because we are concerned that you may be at risk if you use it alone;
- Before hair removal, please shave first if your hair is longer than 1mm. Do not use chemical hair removal method but a shaver. Excessively long hair absorbs heat on the skin surface and can cause skin red and swollen, and hair scorching will produce a peculiar smell, which will affect your experience with the product;
- The equipment is not intended for use by children. Children should be supervised to ensure that they do not play with the equipment to avoid the risk of fire and burns.
- Prevent children from contacting the product to avoid damage, please stay away from children when operating the instrument;
- Do not use on pets to avoid injury;

- Please use different light intensity levels for different body parts as their complexions are different;
- Do not disassemble, repair or modify the device without authorization.
If the device cannot work normally, please send it and related accessories back to the manufacturer for repair by our professional maintenance personnel;
- Power cord can not be used for any other purpose than work. Abnormal use may cause unexpected risks such as choking, electric shock and so on.
- If the power cord is damaged, to avoid danger, it must be replaced by a specialist from the manufacturer's maintenance department or similar department;
- Do not use antiperspirants immediately after irradiation to the underarms;
- Please regularly check and clean the dirt on the light exit to avoid scalding and failure during use;
- Do not place the light exit close to the items that are easy to be demagnetized, so as to avoid malfunction and damage;
- Do not use it when blocking the breather vent, as it may cause smoke and fire;
- Do not point the light exit towards black objects, to avoid the machine running abnormally;
- Do not illuminate clothes, hair, and other items to avoid scorching, discoloration, and damage;
- After use, please wipe off the hair and skin care products attached to the device to avoid malfunction and skin damage;
- Please set the irradiation intensity according to your skin condition. Be sure to use from the first level when using for the first time;
- The product cannot be left unattended after it is powered on, please shut down immediately after use;

- Clean the skin to ensure that there is no hair residue on the skin surface, it is completely dry and free of oily substances and stains; If the skin in the irradiated area is not dry or clean, it will affect the use effect, and may cause the product not to flash or accidentally flash;
- Do not use with essential oil products with the device at the same time. Please ensure that the skin in the irradiated area is clean and dry, and no lotions, lubricants and other cosmetics are used;
- If the ambient temperature difference exceeds 20 degrees (from low temperature to high temperature, or vice versa), please turn off the device and reuse it after about 2 hours;
- Do not use near flammable and explosive materials (including sprayers and gas stations, etc.), which may cause fire and electric shock;
- The ambient temperature for this hair removal device is 0-35°C, please confirm the room temperature before use;
- Within 24 hours after irradiation, please avoid direct ultraviolet rays;
- Do not expose the product to direct sunlight or ultraviolet light for a long time, which will seriously damage the product;
- If the built-in fan dose not work (no murmur is heard) after turning on the power, do not use the device to avoid malfunctions such as smoke and fire;
- Please use this product with caution for those who have used chemical hair removal within 30 days, those who have used skin whitening products recently, or those who have used cosmetics containing vitamin A ingredients;
- Do not use for the following body parts or the following conditions:
 1. Intra-oral, genitals, mucosa, etc.
 2. Body parts that have undergone plastic surgery, surgery
 3. Throat knot

4. Eyeball, eyelid, around eyes

5. Head

6. Acne area

7. Spots and moles

8. Eyebrows, hair

9. Ear

10. Nipple, areola, penis

11. Inside the navel

12. Large vein

13. Tattoo, body painting

(Note: This list is not exhaustive, if you are not sure whether you can use this product, we recommend that you consult a doctor).

- Do not use the device on naturally dark skin;
- Don't expose yourself to excessive sunlight before hair removal;
- When preparing hair removal, make sure the light exit is close to the skin, and then perform the hair removal operation;
- It is prohibited to use the device in dark environment ;
- Do not use the device in direct sunlight to avoid device failure;
- Do not violently impact or shake this product;
- Do not unplug directly without turning off the device to avoid damage to the device;
- When stop using the hair removal device for a break, please make sure that the LCD is facing left or right, which is beneficial to thermal dissipation and cooling.
- Men can not use this device on the face.

1.3 Adverse reaction

The skin may become slightly hot and red after use, which is normal and it will soon disappear. But you should stop using the device immediately if your skin blisters, burns or is allergic, as continue using may cause burns or other skin damage.

2 Product introduction

2.1 Indications for use

IPL Cooling Hair Removal Device is an over-the-counter device intended for removal of unwanted body and/or facial hair.

2.2 About IPL Cooling Hair Removal Device


Innovative sapphire cooling skin care technology can effectively cool the skin, reduce skin heat, avoid skin burn and better protect skin. Sapphire plane head is compressed against the skin which presses down the hair follicles to make its roots closer to the areas with more energy, thus improving the efficiency of energy absorption of hair follicle. Sapphire plane lamp head is easy to clean and has good transparency which brings a better hair removal effect for its stronger light transmittance.

The product sets 6 energy levels to meet the needs of different skin sensitivity. The device and power adapter can be separated, which makes the device compact and portable, easy to store. This product has an integrated lamp holder that does not need to be replaced for a long life, and is suitable for multiple hair removal areas (such as: face, lips, underarms, bikini lines, arms, legs, etc.), which can save expensive lamp holder replacement costs, and you have no need to buy accessories frequently and there are no hidden additional costs. In addition, it has a perfect filtering system, which can filter ineffective short wavelength and harmful light very well.

2.3 Packing list

Upon opening the product package, you will find the following contents inside:

- IPL Cooling Hair Removal Device Host ×1
- Adapter ×1
- Goggles ×1
- Razor×1
- User Manual ×1
- Quick Start Guide ×1
- The Letter Just for U ×1

IPL Cooling Hair
Removal Device Host

Adapter

The Letter Just for U

Goggles

User Manual


Razor


Quick Start Guide

2.4 Product components

LCD screen

3 How it works

3.1 How Does IPL Work

IPL Cooling Hair Removal Device is a Pulsed Light system which emits intense pulsed light (IPL) at a wavelength ranging from 550-1200nm. The device works on the principles of selective photothermolysis. That is, thermal damage will be triggered to target chromophores by using light of appropriate wavelength in pulses that exceeds the chromophores' thermal relaxation time but sparing normal skin by limiting the pulse width below the thermal relaxation time for skin.

3.2 Treatment Schedule

Use 3 times a week for weeks 1-2; use twice a week for weeks 3-5; use every half or one month after 5 weeks.

Note: Please use it multiple times according to the growth of the hair. The specific situation varies from person to person

4 How to use

4.1 Application range

RECOMMENDED TREATMENT AREA

IPL Cooling Hair Removal Device is an over-the-counter device intended for removal of unwanted body and/or facial hair.

USE ON ACCEPTABLE SKIN TYPES

This device is to only be used on Fitzpatrick Skin Types I-V.

BEFORE USING THE DEVICE, choose your skin tone on the basis of Fitzpatrick Phototyping Scale, which was developed by Harvard Medical School dermatologist, Thomas Fitzpatrick, M.D., PhD, to classify a person's complexion and their tolerance of sunlight. There are six levels of skin tone in the scale.

Please IDENTIFY your skin tone and check the following diagram to see whether the product is applicable:

Safety grade	Applicable				Applicable, but Caution	Not Applicable
Skin type /color	Type I Pale white skin	Type II white skin	Type III Cream white skin	Type IV Moderate brown skin	Type V Dark brown skin	Type VI Deeply pigmented dark brown to darkest brown, black
Example						

Cautions:

- Choose the comfortable level according to your skin tolerance and increase level from one to five levels gradually. If there is no discomfort or reaction, you could increase level gradually.
- The color of your skin may change, for example because your skin has been tanned, it may be darker than the skin tone during the previous treatment. Therefore, it is important that you check the settings to be used before each treatment and set the settings according to different treatment areas. Some areas may be darker than others.
- User with skin type V should use the device with a particular caution because you are more likely to experience adverse events. If you have trouble identifying yourself as Skin Type V versus VI, please consult your physician.

USE ON ACCEPTABLE HAIR TONE

Illustration							
Color	White / Gray	Blonde	Red	Light Brown	Brown	Dark Brown	Black
Applicable or not	NO			Yes			

4.2 Patch test

Before actual treatment, please refer to the Recommended SKIN TONE Table (Section 4.1). Conduct a skin test on a small area 24 hours in advance to check whether there will be any adverse reaction on the skin.

HOW TO TEST

- Choose a small area on the arm. Conduct the skin test slowly from low to high energy levels according to skin tolerance (1~2 shots). If there are any signs of swelling, burn mark, change in hair follicle, pain, or other unsuitable symptoms, please do not use this device until a thorough consultation from a physician.
- If there is no skin abnormality after 24 hours, but you feel intense hotness during the treatment, then the last lower level is the most proper energy for your skin.

4.3 Operating Steps

Preparation before use

1. On the basis of Fitzpatrick Phototyping Scale, identify whether your skin tone is suitable for use;
2. Prior to use, make sure your skin is in good condition without any injury and wound. Use a razor to shave off the hair in the area in need of hair removal. Thoroughly clean the skin to ensure that it is hairless, completely dry, and free from oily substances (residues such as perfume, lotion, deodorant and sunscreen).

⚠ The razor should not touch the skin directly. Use it with foam or shave while taking a shower. The razor is sharp. Please use it carefully to avoid scratches.


3. Check that there is no dust, hair, or other dirt remaining on the light exit, and that there is no damage or abnormality. If there is dirt, please clean it to avoid product failure, deterioration and skin problems

Product operating steps

- 1** Use a razor to shave the targeted area for hair removal.


- 2** Plug the power cord into the power input connector on the rear of the product. After connecting the power cord to the product's power input connector, insert the plug into a power outlet.


- 3** Wear the goggles before you conduct hair removal operation to protect your eyes.


- 4** Long press the Power On-Off key for 2s on the back of the host to turn on the device, the first level will be displayed when the device is turned on, in the meantime, the cooling function will be activated and the cold compress indicator will light on.


- 5** Short press the Power On-Off key to adjust the level position. The device has six level positions. Adjust the required hair removal level position. The higher the level, the greater the energy, the faster the effect. It is recommended to start from the first level for the first use, then step up to acceptable level (dark skin is prohibited to use)

- 6** Select the continuous flash mode or single flash mode according to the target area and place the light exit on the skin to be treated and apply it vertically. Then enjoy your journey of hair management!

- 7** After hair removal is completed, press and hold the Power On-Off key for 2s to turn off the device (never pull the plug without turning off the device to avoid damage)

Mode introduction

The only difference between continuous flash mode and single flash mode is that the continuous flash mode controls the device to automatically emit light pulse at a certain interval without manual operation, while under single flash mode, user needs to manually press the flash key. Change the flash modes will not change the technical specifications of the device.

· Continuous flash mode (for large areas such as arms, legs):

Long press the flash key for 2s, the continuous mode is turned on, the level remains unchanged after switching. Place the light exit on the area to be treated on the skin and fully fit the skin, to flash light automatically. Immediately start moving the device slowly and uniformly on the skin, and the flash will be emitted at a certain interval. At this time, please keep the hair removal device moving, and slide back and forth twice in the area where hair removal is needed, to ensure that the hair removal area is covered evenly and avoid missing light. To prevent over-dose, please do not irradiate the same skin area for more than 3 times continuously at a time.

Note: The movement speed can be adjusted according to the flash speed, so that the coverage is more uniform. The flash frequency is faster in the first level, and it is slower in the second level. The higher the level, the slower the light output speed.

· **Single flash mode (suitable for small areas such as fingers and lip hair):**

Long press the flash key under continuous flash mode to switch to single flash mode, the level remains unchanged after switching. Place the light exit on the area to be treated on the skin and fully fit the skin. At this time, the hair removal sensor light will light up, indicating that the lamp head is in full contact with the skin. Short press the flash key to start hair removal. Move the lamp head to the next skin area, wait for the green light of the hair removal sensor light to light up, and then press the flash key to repeat the above operation in order.

Caution:

At any time during the hair removal process, if your skin begins to feel over-heated, just pause for a few seconds, allow the light exit to cooling down again, and then perform hair removal.

Over-temperature protection:

When the internal temperature of the device exceeds 60°C for 30 seconds, the over-temperature protection function of the device will be triggered which will cause the device not to emit the treatment light pulses. When the temperature drops to 50°C for 30 seconds, the device will start working again.

Cooling function:

1. The continuous flash of the hair removal device will continuously accumulate heat. In order to keep the cooling effect effective, it is recommended to use the sapphire lamp head to cool down the irradiated skin area for about 30s after hair removal, and then perform the hair removal operation in the next area. This will cool down the skin again and make sapphire cooling effect better.

Note: The cooling function has no medical therapeutic effect and is just for cool down the skin temperature.

2. Press the lamp head tightly against the skin to make sure that the lamp head is in full contact with the root of the hair follicle and the hair follicle will absorb higher energy, then achieve a better hair removal effect. Sapphire cooling skin care technology can effectively cool the skin, reduce skin heat and better protect skin.

4.4 Skincare after treatment

Cautions after hair removal

- It is recommended to avoid contact with water within 6 hours after the hair removal operation, and then moisturizing lotion can be used to moisturize the irradiated area.
- Do not use perfumes, facial mask, or cosmetics containing buffing beads after skin hair removal
- Avoid UV exposure for 24 hours after hair removal, and apply sunscreen if necessary.
- Sunbathing is likely to cause skin darkening or sunburn. Irradiation will cause redness, swelling, and inflammation. Therefore, during the use of this product, be careful not to be exposed to the sun, since the time to repair the sunburn is not less than 1 week.
- Artificial tanning is not allowed during the hair removal operation. If you need plastic surgery after hair removal (only for the hair removal area), consult your doctor.
- If you feel unwell, consult your doctor.

5 Troubleshooting

This chapter outlines the problems you may encounter when using this product. If you cannot solve the problems encountered with the following solutions, please contact the customer service center in your country/region.

- When an abnormal situation occurs or the fault cannot be eliminated, please contact the after-sales customer service in the corresponding language.
- Please keep your purchase invoice and warranty card, you will need to bring them when your product fails

Problem	Possible reason	Solution
Unable to boot normally	Power is not connected	Check if the power is plugged in
Unable to emit light	Whether the temperature of the sapphire lamp head and the body is too high; Whether the head is close to the skin	Wait for the device to cool down before using it. Hold the host close to the skin until the light indicator lights up
The work indicator is off	Bad device contact	Check if the power supply is working properly; Repair or return product immediately
Odor appears when flashing	Filter is contaminated; Hair residue on body surface	Wipe filter; Shaving hair and reducing light intensity

Unbearable pain during hair removal	Did not shave off the hair on the hair removal area	Before using this product, please shave off the hair on the area to be irradiated. If shaving causes skin irritation, cut your hair as short as possible and use it in a comfortable setting.
	Hair removal site is prohibited to use this product	Do not use this product for inner labia, vagina, anus, nipples and areola. Men should not use this product on the face or penis
	Light intensity is too high	Reduce light intensity to a comfortable level
	The light exit is damaged	Repair or return product immediately
After using, the hair removal area turns red	A slight reddening is normal to the body, and it will soon fade	Use aloe vera gel to soothe the skin or use a sapphire lamp head for cooling
Unsatisfactory hair removal effect	You use a light intensity level that is too low for you	Please choose a higher level next time
	Missing flash areas on the hair removal area	Flashes in close proximity and there should be some overlap between the flash areas
	This device has no effect on your hair color	If your hair is pale blonde, gray, red and white, the exposure will not be effective
	Insufficient frequency of use of this device	To successfully remove all hair, we recommend that you follow the recommended irradiation regimen
	Your response to light exposure is slow	This product should be used continuously for at least 5 months, as it can still suppress hair regeneration during this period
The product becomes hot and humid during use	Maybe due to long-term device operation	Turn off the device and let it cool down before using it

6 Technical parameters and specifications

Basic specifications			
Product name	IPL Cooling Hair Removal Device		
Model	HR-544		
Dimension	9.3*3.0*1.7 inches		
Weight	0.74 lbs		
Type of Protection	BF Type, Class II		
Software version	V1		
Output specifications (The parameters are the same under single/ continuous flash mode)			
Type of light	Intense Pulsed Light (IPL)		
Light wavelength	550-1200nm		
Pulse duration	6.4-7.2ms		
Interval time between flashes	1.3s~2.6s		
Treatment area	3.6 cm ²		
Treatment level	Level 1-6		
Energy density (±20%)	2-4 J/cm ²		
Output energy (±20%)	Level 1: 7.14J Level 4: 12.48J	Level 2: 10.0J Level 5: 12.88J	Level 3: 11.86J Level 6: 13.62J
Power specifications			
Power supply	Supplied by external power adapter Input: 100-240V, 50/60Hz Output: DC12V 4A		
Power	48W		
Usage Specifications			
Expected service life	More than 500,000 counts		
Operating conditions	Operating temperature: 0-35°C Humidity: 10%-90% RH Barometric pressure: 86-106kPa		
Storage and transportation conditions	Temperature: -10-60°C Humidity: 5%-90% RH Barometric pressure: 50-106kPa		

7 Maintenance

Cleaning

In order to ensure the best results and longest service life of the product, please clean the product before and after each hair removal to avoid long-term accumulation of hair and dust. If you do not clean it properly, the product will not get its proper effect.

- Before starting cleaning, the device must be switched off and unplugged to allow the device to cool down.
- Use a soft, clean and dry cloth dipped in water to wipe and clean the mirror surface of the light outlet.
- Never wash the product and accessories in faucets and dishwashers.
- Do not use steel wool, abrasive cleaners (toothpaste, etc.) or corrosive liquids (gasoline, motor oil, acetone, etc.) to clean the product.
- If water does not clean well, you can use cloth dipped in 75% medicinal alcohol to clean.

Storage

- To avoid breakage or missing device and its accessories, please store IPL Cooling Hair Removal Device and its accessories in the gift box after each use.
- Clean the device before you store it.
- Store the device in a dry, dust-free environment at -10~60°C.

8 Disposal

In order to avoid environmental pollution, please dispose this device according to local environmental requirements and do not discard casually.

9 EMC statement

1. IPL Cooling Hair Removal Device needs special precautions regarding EMC and need to be installed and put into service according to the EMC information provided in the accompanying document;

2. Portable and mobile RF communications equipment can affect IPL Cooling Hair Removal Device.

Warning:

1. Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

2. Use of accessories, transducers and cables other than those specified or provided by the manufacturer of IPL Cooling Hair Removal Device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

3. Use of IPL Cooling Hair Removal Device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

10 Symbol meaning

Symbol	Meaning	Symbol	Meaning
	Batch Code		Serial number
	Refer to instruction manual/ booklet		Type BF Applied Part
	Class II equipment		Eye protection
	Manufacturer		Date of manufacture
	Caution		General warning sign
	Keep away from sunlight		Keep dry
	Fragile, handle with care		“WEEE (Waste Electrical and Electronic Equipment)” .The waste products should be handled legally.

11 Warranty

- The warranty period of this product may vary according to the country in which it is sold, and we will try to extend it as far as possible under the laws of that country.
- The warranty period for this product is 1 year from the date of purchase. The power supply and accessories are guaranteed for one year from the date of purchase.
- If the device has quality problems during the warranty period, Jizhimei will repair or replace it free of charge. However, device damage caused by improper use or improper operation is not covered by the warranty. Jizhimei will replace the damaged part or the old one at the cost price. The old device is owned by Jizhimei.

- Please present your warranty card when you need a warranty. The Jizhimei product's warranty date is based on the date when you purchased it. If you do not have a purchase date, it will be calculated based on the factory date.

In the following cases, Jizhimei may not provide warranty or return.

Jizhimei is not responsible for reissuing the purchase invoice and warranty card for you, please keep your purchase invoice and warranty card for yourself.

1. For product damage caused by unpredictable events such as fire, earthquake, and flood.
2. Device failure and damage caused by your negligence, misuse, or use under other abnormal conditions.
3. Damage to the device caused by improper disassembly and storage (such as rodent damage or liquid penetration).
4. Damage caused by mandatory use of the product in a non-product use environment, including but not limited to exceeding the workload.
5. Product damage caused by other conditions that Jizhimei cannot control.
6. For the function of this product or the incidental damage to its use (such as loss of your business benefits, business interruption, other losses caused by the use of the device, the cost of using this product to replace equipment, devices or services, your customer claims etc.), Jizhimei assumes no responsibility for this.
7. For the type and extent of the accident, Jizhimei will bear it according to relevant laws and regulations.

12 Contact information

NAISIGOO × | 📞 🌐

WhatsApp Tel: +86 18938936202

Shenzhen Jizhimei Technology Co., Ltd.

Address: Building A203, Building A/B, Zone B, Gonghe Community,
Shajing Street, Baoan District, Shenzhen

Appendix A. EMC information

Essential performance:

The energy density of IPL Cooling Hair Removal Device should be in the range of 2-4 J/cm², error is ±20%.

Table 1

Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Group B
Harmonic emissions IEC 61000-3-2	Group A
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies

Table 2

Guidance and Manufacturer' s Declaration – Electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact, ±2 kV, ±4kV, ±8 kV, ±15 kV air	
Electrical fast transient /burst IEC 61000-4-4	±2 kV for a.c. power cord	
Surge IEC61000-4-5	±0.5kV and ±1kV for a.c. power cord to power cord	
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT; 0.5 cycle at 0°, 45°, 90°,135°, 180°, 225°, 270°, 315° 0 % UT; 1 cycle 70 % UT; 25/30 cycle 0% UT; 250/300 cycle	
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	
Conducted RF IEC61000-4-6	3 V rms. 0.15 MHz to 80 MHz 6 V rms in the ISM and amateur bands between 0.15 MHz and 80 MHz	
Radiated RF	10 V/m	
IEC 61000-4-3	80 MHz – 2,7 GHz 80 % AM at 1 kHz	
NOTE UT is the a.c. mains voltage prior to application of the test level.		

Table 3

Guidance and Manufacturer' s Declaration - IMMUNITY to proximity fields from RF wireless communications equipment					
Immunity Test	IEC60601 test level				Compliance level
	Test frequency	Modulation	Maximum power	Immunity level	
Radiated RF IEC 61000-4-3	385 MHz	**Pulse Modulation: 18Hz	1.8W	27 V/m	
	450 MHz	*FM+ 5Hz deviation: 1kHz sine	2 W	28 V/m	
	710 MHz 745 MHz 780 MHz	**Pulse Modulation: 217Hz	0.2 W	9 V/m	
	810 MHz 870 MHz 930 MHz	**Pulse Modulation: 18Hz	2 W	28 V/m	
	1720 MHz 1845 MHz 1970 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	
	2450 MHz	**Pulse Modulation: 217Hz	2 W	28 V/m	
	5240 MHz 5500 MHz 5785 MHz	**Pulse Modulation: 217Hz	0.2W	9 V/m	
Note* - As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be the worst case. Note** - The carrier shall be modulated using a 50 % duty cycle square wave signal.					

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